

PerfeCta

repliQa

Q

qScript

sparQ

Extracta

AccuStart

Sample
Preparation

Reverse
Transcription

PCR & qPCR

Genotyping

microRNA
Profiling

NGS

CONTENT

1.0	Quantabio Brand Story	4
1.1	About Quantabio	4
1.2	Our History	5
1.3	Markets & Applications	6
1.4	Family of Brands	8
1.5	Core Technologies	9
2.0	Sample Preparation	12
2.1	Extracta DNA Prep for PCR	12
2.2	Extracta DBS	13
2.3	5PRIME Phase Lock Gel	15
3.0	Reverse Transcription	17
3.1	First-strand cDNA Synthesis	17
3.1.01	qScript® cDNA SuperMix/qScript XLT cDNA SuperMix	18
3.1.02	qScript cDNA Synthesis Kit	20
3.1.03	qScript Flex cDNA Synthesis Kit	20
3.1.04	PerfeCта DNase I	20
3.2	Reverse Transcription PCR (RT-PCR)	21
	Conventional RT-PCR	
3.2.01	qScript XLT One-Step RT-PCR Kit	22
	Quantitative RT-qPCR	
3.2.02	Qscript Iyo 1-step	24
3.2.03	UltraPlex One-Step ToughMix	26
3.2.04	qScript One-Step SYBR® Green RT-qPCR	27
3.2.05	qScript XLT One-Step RT-qPCR ToughMix®	29
3.2.06	qScript One-Step RT-qPCR	31
4.0	PCR & qPCR	32
4.1	qPCR Instrumentation	32
4.1.01	Q qPCR Cyclers	32
4.2	Conventional PCR	34
4.2.01	AccuStart™ Long Range SuperMix	34
4.2.02	repliQa™ HiFi ToughMix	36
4.2.03	AccuStart II PCR ToughMix	40
4.2.04	AccuStart Taq DNA Polymerase HiFi	42
4.2.05	AccuStart II GelTrack® PCR SuperMix	43
4.2.06	AccuStart II PCR SuperMix	44

4.2.07	AccuStart II Taq DNA Polymerase	45
4.2.08	5PRIME HotMaster Taq DNA Polymerase	46
4.2.09	5PRIME HotMasterMix	47
4.2.10	10 mM dNTP Mix	48
4.3	Real-Time Quantitative PCR	49
	Multiplexed Pre-Amplification		
4.3.01	PerfeCta® PreAmp 5x SuperMix	50
	SYBR® Green Detection		
4.3.02	PerfeCta SYBR Green FastMix®	52
4.3.03	PerfeCta SYBR Green SuperMix	52
	Probe-based Detection		
4.3.04	PerfeCta qPCR ToughMix	55
4.3.05	PerfeCta FastMix II	58
4.3.06	PerfeCta MultiPlex qPCR ToughMix	59
4.3.07	PerfeCta MultiPlex qPCR SuperMix	61
4.4	Cloning	62
4.4.02	repliQa HiFi ToughMix	62
4.4.02	repliQa HiFi Assembly Mix	66
5.0	Genotyping	68
5.1.01	AccuStart Genotyping ToughMix	68
5.1.02	AccuStart II PCR Genotyping Kit	70
5.1.03	AccuMelt HRM SuperMix	71
6.0	microRNA Profiling	73
6.1.01	qScript microRNA cDNA Synthesis Kit	73
7.0	Next Generation Sequencing (NGS)	77
7.1.01	sparQ DNA Frag & Library Prep Kit	77
7.1.02	sparQ DNA Library Prep Kit	81
7.1.03	sparQ HiFi PCR Master Mix	85
7.1.04	sparQ PureMag Beads	87
7.1.05	sparQ Fast Library Quant Kit	90
7.1.06	PerfeCta NGS Quantification Kit	94
	Instrument Compatibility	96
	License Disclaimer and Trademarks	97
	Contact	99

1.1 ABOUT QUANTABIO

A fast-growing, innovative brand and leading provider of advanced DNA and RNA amplification reagents for the most demanding molecular testing applications in applied, translational and life science research.


Our Recipe for Success

It's in our DNA. Be Resilient. **ADAPT.**

It starts from ideation with product design & development, continuing all the way to customer application support. We aim to make and **ADAPT** product solutions that deliver superior results, enable easier workflows and better affordability while providing unwavering support in ever changing life science and molecular diagnostic environments.

Your success is our success.

How we **ADAPT** product solutions in 4 Steps



Step 1

Start with
Enzymes
Engineered for
specific application


+



Step 2

Formulate with
Tough Additive
To optimize reagent
compositions

=



Step 3

Create unique
Superior
Performing
products

Applications

- RT
- qPCR
- PCR
- NGS

Step 4

Provide unwavering
Customer Support
For multiple molecular
biology workflows




1.2 OUR HISTORY

Legacy of Innovation

The Quantabio team leverages decades of experience in developing pioneering amplification technologies to deliver cutting-edge products to researchers focused on critical cloning, PCR, qPCR and Next-Generation Sequencing (NGS) based applications.


Founded 2001 by former amplification scientists from Life Technologies (Invitrogen)

Quanta Biosciences founded



1.3 MARKETS & APPLICATIONS

Quantabio utilizes its novel DNA and RNA amplification technologies along with proprietary buffer chemistry to create unique and differentiated reagent solutions for life science applications ranging from sample extraction to cDNA synthesis, genotyping, gene expression, cloning and even next generation sequencing. Our products are used in thousands of labs around the world for life science research as well as in applied testing markets.




Translational Research Markets

Quantabio reagents are the gold standard for sensitive and reliable quantitative assay performance in PCR, qPCR and NGS. As an example, our proprietary ToughMix® reagents enable efficient amplification for challenging PCR inhibitors.

Quantabio products are used globally for detection of a variety of infectious diseases and in new born screening assays such as:

- Polio
- Influenza
- West Nile
- Zika
- Spinal Muscular Atrophy (SMA)
- Severe Combined Immunodeficiency (SCID)

Genetic Mutation Detection



Pathogen Detection



Applied Markets

Quantabio products are used routinely for applications such as:

- GMO Testing
- Bacterial Contamination
- Beer Spoilage

Food Testing



Animal Health



Environmental Testing



Plant Testing



1.4 FAMILY OF BRANDS

Extracta

Extracta extraction reagents provide a simplified and cost-effective alternative to traditional nucleic acid (NA) purification methods and are optimized to work in series with ToughMix products. Optimized for clinical specimens (blood and dried blood spots), plant and animal tissues and environmental samples.

AccuStart

AccuStart ultrapure DNA polymerase, contains a stringent antibody hotstart to ensure specific and efficient primer extension only after activation at 94°C and rigorously purified to remove host *E. coli* genomic DNA as residual host DNA in typical recombinant enzyme preparations can limit assay sensitivity and obscure detection of low copy targets.

PerfeCta

PerfeCta qPCR reagents combine a stringent, ultrapure antibody hotstart with performance engineered DNA polymerase in stabilized 1-tube formulations optimized for the specific performance needs of real-time quantitative PCR. Proprietary additives help eliminate persistent bubbles to enable efficient vortex mixing and fewer technical replicates thereby conserving precious sample. Adaptive buffer chemistry accommodates most assay designs and can be used with existing assay designs.

qScript

qScript reverse transcription reagents leverage proprietary, performance-engineered qScript® reverse transcriptase in a variety of stabilized, user-friendly reagent formulations that maximize cDNA yield and provide linear cDNA synthesis across a broad dynamic range of input RNA. qScript cDNA synthesis kits redefine what is possible in speed, convenience, reproducibility, specificity and limit of detection (LOD) sensitivity in qPCR and RT-PCR applications.

repliQa

repliQa HiFi Assembly Mix & repliQa Hi Fi ToughMix simplifies the construction of recombinant DNA through seamless assembly of multiple DNA fragments in a single, isothermal reaction. Similar in principle to the Gibson Assembly® Method, the high efficiency repliQa mixes are ideal for a range of genetic engineering applications: cloning, site-directed mutagenesis and synthetic biology.

sparQ

sparQ offers complete solutions for library preparation, amplification, purification and quantification for Illumina next generation sequencing platforms. High quality reagents deliver unmatched efficiency and robust performance to ensure reliable and reproducible sequencing results while reducing total sequencing costs.

1.5 CORE TECHNOLOGIES

Manufacturing Excellence

- ISO 13485 quality certified
- Ultrapure, performance-engineered enzymes
- Ultralow residual host *E. coli* DNA

Engineered Stability

- Stringent enzyme activation control with AccuStart antibody technology
- Reaction setup and multi-day storage at ambient temperature
- Impervious to repetitive freeze-thaw

Formulated for Quantitative Real-Time PCR

- Optimized 1-tube reagents minimize pipetting steps and improve accuracy
- Supports efficient vortex mixing and eliminates error-causing persistent bubbles
- Inert AccuVue plate loading dye provides visual confirmation of reaction assembly

Tough-Tested

- ToughMix reagents withstand a broad spectrum of PCR inhibitors
- Reliable assay performance with challenging starting materials and crude extracts

Optimized to improve sequencing performance & economics

- Comprehensive solutions for DNA fragmentation, library preparation, amplification and quantification
- Novel formulations streamline NGS workflows reducing total turn-around-time
- Proprietary enzymes & buffer compositions improve library yields and sequencing results from low inputs

RT-qPCR is a powerful molecular biology technologies.



RT-qPCR

Reproducible

Quantabio kits define experimental reproducibility through multiple proprietary technologies: ultra low *E. coli* DNA, low foam, and enhanced stability. Due to our exacting lot-to-lot consistency standards and innovative 1-tube formulations that minimize pipetting, our reagent technologies provide highly consistent results. Patented additives actively reduce intra-assay variability and withstand repetitive cycles of freeze-thaw to deliver assay reliability.



RT-qPCR

ToughMix

PCR inhibitors are common in crude samples and readily compromise assay performance. Use of Quantabio TOUGH-tested ToughMix reagents results in enhances PCR performance with crude or contaminated samples. The advanced ToughMix buffer technology is engineered to withstand several types of PCR-inhibition, providing robust and reliable results from a variety of starting materials and purification methods.



RT-qPCR

qScript

We put the "Q" in RT-qPCR with our advanced reverse transcriptase technology that is synonymous with maximum yield and sensitivity. qScript first-strand cDNA synthesis reagents are rigorously optimized to provide sensitive and reliable detection of low abundance RNA for qPCR assays. The broad, linear dynamic range of input RNA (10 pg - 1 µg) provides reliable assay sensitivity for robust gene expression analysis.

tool that represents the Quantabio core



RT-qPCR

Perfecting qPCR

Our PerfeCTa real-time quantitative PCR reagents are rigorously optimized; all-in-one reagents that dramatically simplify reaction setup and contain patented technologies to actively reduce assay variance. Our robust, ultrapure antibody hot start AccuStart technology drives precise target amplification with the absolute maximum limit of detection sensitivity.



RT-qPCR

Customization

Our proprietary formulation processes allow us to rapidly configure customized reagent solutions according to client specific needs. Whether it's a defined lot size for a large scale project, a packaging or fill volume requirement to suit a particular workflow, or modified composition to fine tune assay performance, Quantabio has the flexible responsiveness to realize your custom needs.



RT-qPCR

Reliable

Quantabio is your trusted reagent supply partner. We pride ourselves on manufacturing quality excellence and an industry-leading technology portfolio.

2.0

Sample Preparation

Extracta DNA Prep

Quick and easy DNA extraction of PCR-ready genomic DNA

Extracta DNA Prep is an entirely reagent-based system for extracting and stabilizing template DNA from a variety of biological starting materials for sensitive downstream applications such as PCR, qPCR and HRM analysis


FEATURES AND BENEFITS:

- Simple, reagent-based system requires minimal technical skill
- Incubation step can be carried out in 96-well PCR plates or tubes using a standard DNA thermal cycler
- Compatible with a wide-range of clinical specimens, plant and animal tissues, and environmental samples
- Optional stabilization buffer allows for extended storage of extracted DNA templates

DESCRIPTION:

Extracta DNA Prep for PCR is a two-component reagent kit for rapid extraction of PCR-ready genomic DNA from a variety of tissues. Samples are processed in less than 30 minutes with minimal hands-on time and technical skill. Extracted genomic DNA is suitable for sensitive downstream PCR applications including end-point PCR, High Resolution Melt Analysis (HRM) and quantitative real-time PCR (qPCR) without requiring any additional clean-up. In addition, the extracted DNA may be used in multiplexed PCR applications such as transgene or knock-out analyses. Tissue extractions can be done in tubes, plates or deep-well blocks to allow for adaptation to workflow and automation on liquid-handling workstations.

Extracta Procedure



ORDER INFO

Product Name

Extracta DNA Prep - 2.5 ml
Extracta DNA Prep - 25 ml
Extracta DNA Prep - 250 ml

Quantabio Catalog Number

95091-002
95091-025
95091-250

Size

2.5 ml
25 ml
250 ml



Extracta DBS

PCR-ready genomic DNA from dried blood spots

FEATURES AND BENEFITS:


- Optimized for DNA extraction from dried blood spot punches
- Single reagent for PCR-ready DNA in 30 minutes
- Maximized assay sensitivity, lower Cq values, when combined with Quantabio ToughMix
- Compatible with high-throughput automation for PCR, qPCR and NGS applications

DESCRIPTION:

Extracta DBS is a ready-to-use DNA extraction reagent for rapid and efficient recovery of PCR-ready DNA from dried blood spots (DBS) on Guthrie cards or Whatman 903 filter paper. This patented single-solution process produces DNA eluates that are substantially free of PCR inhibitors and compatible with a variety of end-point PCR, real-time PCR and Next Generation

Sequencing (NGS) or Sanger Sequencing reagents. Application of Extracta DBS with PerfeCta qPCR ToughMix or PerfeCta MultiPlex qPCR ToughMix enables accurate and reproducible quantification of DNA sequences in blood using TaqMan hydrolysis probe real-time qPCR.

Fast, single solution extraction workflow in 30 min




2.1 Extracta DBS DNA extraction from dried blood spots workflow. PCR-ready DNA is eluted from dried blood spots with less handling in 30 minutes.



Combine with ToughMix and maximize qPCR assay sensitivity

Extracta DBS is the perfect match with Quantabio ToughMix for sensitive and precise target quantification. The crude extraction combined with ToughMix, a Quantabio master mix that is tolerant to common PCR inhibitors, results in higher DNA yields independent of DNA sample inputs and qualities to enable accurate detection and high sensitivity even with low copy targets.

Extracta DBS TREC Quantification




2.2 Illustrates the results of a T-cell Receptor Excision Circles (TREC) assay using Extracta DBS and PerfeCra ToughMix. Samples were generated using dried blood spot punches following the Extracta DBS protocol and used subsequently for quantification of T-cell Receptor Excision Circles. The samples are representatives of High, Medium and Low TREC copy numbers along with the corresponding Cq values.

Reliable, Consistent Lot-to-Lot Performance

Low amounts and quality of DNA recovered from dried blood spots commonly restrict the utilization of DNA. To overcome these limitations, Extracta DBS increases the yield and quality allowing for efficient and reliable recovery of DNA from dried blood spots. Manufactured under stringent ISO 13485 standards, Extracta DBS ensures uniform lot-to-lot performance resulting in reliable reproducibility in combination with Quantabio ToughMix.

Extracta DBS TREC quantification: uniform lot performance



2.3 This figure demonstrates consistent lot-to-lot performance in a TREC quantification assay using genomic DNA extracted from dried blood spots. Lot-to-lot performance was tested for High, Medium and Low copy number TREC samples. The results highlight the reliability and reproducibility across various product lots which are attributed to Quantabio's high manufacturing and production standards under ISO 13485.

ORDER INFO

Product Name

Extracta DBS, 10 ml
Extracta DBS, 500 ml

Quantabio Catalog Number

95171-010
95171-500

Size

10 ml
500 ml

Quantabio ToughMixes are also available with different concentrations of ROX and in larger reaction sizes.



5PRIME Phase Lock Gel

Simplifies organic extraction of nucleic acid template and improves safety

Organic extraction methods are cost-effective and result in the highest yields of nucleic acid template but are not user friendly and involve hazardous chemicals

FEATURES AND BENEFITS:

- Eliminates interphase contamination of nucleic acid solution
- 30% greater yield of nucleic acids over conventional method
- Gel barrier allows easy sample decanting
- Reduced contact with hazardous organic solvents

DESCRIPTION:


Phase Lock Gel (PLG) is a unique product that eliminates interphase-protein contamination during phenol extraction of DNA or RNA. PLG reduces hands-on time and improves nucleic acid recovery. PLG migrates under centrifugal force to form a tight seal between the aqueous and the organic phase. The organic phase and the interphase materials are effectively trapped in or below the barrier. The stable barrier enables a complete and easy transfer of the aqueous, nucleic acid containing upper phase to a fresh tube.

The benefits are increased yields by up to 30%, better protection from exposure to hazardous compounds and no risk of sample contamination with interphase debris. PLG can be adapted to virtually any protocol requiring extraction of an aqueous sample with phenol and/or chloroform. For convenience, PLG is provided aliquoted into standard centrifuge tubes of various sizes.

Phase Lock Gel Light and Heavy Applications and Compatibilities

For optimal phase separation the composition of the aqueous phase, the organic phase and PLG must be compatible as the ability of PLG to separate the phases depends on matching

the differences in density between aqueous and organic media. Besides general differences of the organic phase due to the starting material, density is also influenced by salt and protein concentration in the aqueous phase. To ensure compatibility, PLG comes in 2 density formulations, Heavy (H) and Light (L). Choose the formulation that fits your specific application from the table next page:





Phase Lock Gel, Density Selection Chart

Aqueous Phase	Organic Phase			
	PCI	CI	H ₂ O or Buffer saturated PC	H ₂ O or Buffer saturated Phenol
< 0.5 M NaCl	L, H	L, H	L, H	L
< 1 mg/ml BSA	L, H	L, H	L, H	L
Cleared bacterial lysate	H	H	H	–
Plasmid DNA homogenates	H	H	H	–
Tissue homogenates	L, H	L, H	L, H	L
Genomic DNA isolation	L, H	L, H	L, H	L
RNA isolation	H	H	H	–

PCI = 25:24:1 Phenol : Chloroform : Isoamyl Alcohol

PC = 1:1 Phenol : Chloroform

CI = 24:1 Chloroform : Isoamyl Alcohol

– = Conditions unsuitable for Phase Lock Gels

ORDER INFO

Product Name

Quantabio Catalog Number

Size

Phase Lock Gel Heavy 2 ml - 200 Tubes

2302830

200 tubes

Phase Lock Gel Light 2 ml - 200 Tubes

2302820

200 tubes

Quantabio qScript reverse transcriptase technologies set the standard for simplicity, reproducibility, and broad linear dynamic range for quantitative and conventional RT-PCR applications.

qScript cDNA synthesis reagents provide highly sensitive first-strand cDNA synthesis in a variety of easy-to-use reagent configurations for RT-PCR and RT-qPCR.

Ultrapure, performance-engineered M-MLV reverse transcriptases are pre-blended with ribonuclease (RNase) inhibitor protein in rigorously optimized 1-tube SuperMix formulations and separate component kits to suit specific assay designs and workflow preferences.

First-Strand cDNA Synthesis

3.1

PRODUCT OVERVIEW

	qScript cDNA SuperMix	qScript XLT cDNA SuperMix	qScript cDNA Synthesis Kit	qScript Flex cDNA Synthesis Kit
Kit Format	Single Tube	Single Tube	Two Tubes	Five Tubes
RT Enzyme	MMLV, RNase H+	MMLV, reduced RNase H activity	MMLV, RNase H+	MMLV, RNase H+
Priming Method	Oligo(dT) & random primers	Oligo(dT) & random primers	Oligo(dT) & random primers	Oligo(dT), random primers or gene specific primer
RNA Input Range	10 pg – 1 µg	1 pg – 2 µg	10 pg – 1 µg	10 pg – 1 µg
Amplicon Length	1 kb or less	1 kb or less	1 kb or less	12 kb or less
Optimal Reaction Time	40 min	30 – 70 min	40 min	60 – 90 min
Yield (Sensitivity)	+++	++++	++	++



qScript cDNA SuperMix/qScript XLT cDNA SuperMix

Superior cDNA synthesis in a single step

FEATURES AND BENEFITS:

- 5x concentrated SuperMix maximizes input volume with dilute samples of total RNA
- Broad linear dynamic range with limiting (10 pg) or plentiful samples of total RNA
- Pre-blended with ribonuclease inhibitor protein (RIP) and an optimized blend of random hexamer and oligo(dT) primers to ensure accurate representation of 5' and 3' sequences

qScript cDNA SuperMix

Stabilized 1-tube SuperMixes simplify reaction assembly and minimize risk of pipetting error

1x

MMLV RT Buffer, dNTPs, Mg²⁺ Randomers + Oligo-d(T) RNase inhibitor

Add 4 µl SuperMix to RNA template

Incubate

- 1 Tube
- 1 Pipetting Step
- 40 Minutes (70 min for qScript XLT Supermix)

3.1 qScript cDNA SuperMix (left) includes all necessary components in a single tube – just add RNA and go! Compare this to competitor options (right) that introduce numerous opportunities for error.

Other cDNA Kits

Primer

RT Buffer Inhibitor Mg²⁺

Make cocktail on ice-bath

Add primer & lysate to separate tube incubate & snap-chill on ice



Add cocktail to primer/template tube

Incubate

8x




qScript cDNA SuperMix Delivers Higher Yields, Improved Representation of Low Abundance Genes and Superior Linear Dynamic Range



	Slope	Corr.	PCR Eff.	Sensitivity
qScript cDNA SuperMix	-3.625	0.998	88.7%	50 fg
Competitor S	-3.941	0.999	79.4%	500 fg

3.2 Increasing amounts of HeLa total RNA (1 pg – 1 µg) were reverse transcribed using qScript cDNA SuperMix in parallel with another leading supplier kit, according to supplied protocols.

Mix of random primers and oligo dT prevents bias to the 5' or 3' ends of the target





qScript cDNA Synthesis Kit

Economical 2-component kit ideally suited for high throughput gene-expression studies

FEATURES AND BENEFITS:

- Sensitive first-strand cDNA synthesis of RNA sequences ≤ 1 kb for quantitative and conventional two-step RT-PCR
- Broad linear dynamic range with limiting (10 pg) or plentiful samples of total RNA
- Pre-blended with ribonuclease inhibitor protein (RIP) and an optimized blend of random hexamer and oligo(dT) primers to ensure accurate representation of 5' and 3' sequences

qScript Flex cDNA Synthesis Kit

Highly flexible first-strand synthesis system suitable for large target RNA sequences

FEATURES AND BENEFITS:

- User choice of RT priming method; oligo(dT), random hexamers, or GSP
- Highly sensitive first-strand cDNA synthesis of large RNA sequences for quantitative and conventional two-step RT-PCR
- Broad linear dynamic range with limiting (10 pg) or plentiful total RNA samples
- Maximize cDNA yield with proprietary Priming Enhancer additive

PerfeCta DNase I

Simple and rapid removal of residual genomic DNA

FEATURES AND BENEFITS:

- Permanent DNase I inactivation for confident first-strand cDNA synthesis with qScript first-strand cDNA synthesis kits
- Single tube



ORDER INFO

Product Name	Quantabio Catalog Number	Size
qScript cDNA SuperMix - 25 R	95048-025	25 rxns
qScript cDNA SuperMix - 100 R	95048-100	100 rxns
qScript cDNA SuperMix - 500 R	95048-500	500 rxns
qScript XLT cDNA SuperMix - 25 R	95161-025	25 rxns
qScript XLT cDNA SuperMix - 100 R	95161-100	100 rxns
qScript XLT cDNA SuperMix - 500 R	95161-500	500 rxns
qScript cDNA Synthesis Kit - 25 R	95047-025	25 rxns
qScript cDNA Synthesis Kit - 100 R	95047-100	100 rxns
qScript cDNA Synthesis Kit - 500 R	95047-500	500 rxns
qScript Flex cDNA Kit - 25 R	95049-025	25 rxns
qScript Flex cDNA Kit - 100 R	95049-100	100 rxns
PerfeCta DNase I - 100 U	95150-100	100 U (1 x 50 µl)
PerfeCta DNase I - 1000 U	95150-01K	1000 U (1 x 500 µl)

Reverse Transcription PCR (RT-PCR)

3.2

PRODUCT OVERVIEW

	Conventional One-Step RT-PCR	Tough-Tested MultiPlex One-Step RT-qPCR	Tough-Tested One-Step Probe-based RT-qPCR	One-Step Probe-based RT-qPCR	One-Step SYBR-based RT-qPCR
Kit	qScript XLT One-Step RT-PCR Kit	UltraPlex One-Step ToughMix	qScript XLT One-Step RT-qPCR ToughMix	qScript One-Step RT-qPCR Kit	qScript One-Step SYBR Green RT-qRT PCR
Detection Chemistry	N/A	Hydrolysis Probes	Hydrolysis Probes	Hydrolysis Probes	SYBR Green I dye
Sensitivity	+++	++++	+++	++	++
Multiplex Compatibility	N/A	>4	<4	<3	No
Reagent Components	2	1	1	2	2
RNA Input (Linear Range)	1 pg – 1 µg	• 1 pg to 100 ng total RNA; • 10 fg to 10 ng poly A(+) RNA; • 10 to 1x10 ⁸ copies viral RNA			
Amplicon Length	4+ kb	<1 kb	<1 kb	<1 kb	<1 kb



Conventional RT-PCR

qScript XLT One-Step RT-PCR Kit

Tough-tested One-Step Reverse Transcriptase PCR (RT-PCR) in a simplified, 2-component reagent system

The qScript XLT One-Step RT-PCR kit provides highly sensitive detection of large, complex RNA in challenging starting materials supporting high-fidelity downstream applications

FEATURES AND BENEFITS:

- ToughMix® reagent technology withstands PCR inhibitors commonly found in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Optional GelTrack dye streamlines workflow for gel electrophoresis
- Temperature stabilized for support reaction assembly at convenient ambient room temperatures
- Preblended with ribonuclease inhibitor protein to preserve RNA integrity during incubation
- 3'-exonuclease proof-reading polymerase supports high-fidelity downstream applications
- Suitable for TA subcloning large RNA sequences exceeding 4 kb in length

DESCRIPTION:

The qScript XLT One-Step RT-PCR Kit is a convenient and highly sensitive 2-reagent system for amplification of complex RNA templates exceeding 4 kb in length. Both enzyme incubation sequences are carried out in the same reaction mixture without opening between procedures. Advanced qScript XLT reverse transcriptase mutant possesses elevated temperature stability and improved template binding affinity for large complex

RNA templates. Supports a wide-range of challenging starting materials with PCR-inhibitor neutralizing ToughMix additives and delivers consistent, reliable assay performance. Ultrapure AccuStart II hot-start Taq DNA polymerase with 3'-exonuclease proof-reading activity provides stringent activation control for sensitive and precise target amplification.



Quantitative RT-qPCR

Qscript Iyo 1-step

Dry. Stable. Easy. Better. 1-step RT-qPCR

FEATURES AND BENEFITS:



Lyophilized Single Tube Format – Easy-to-use, reduces cross-contamination



High Sensitivity & Specificity – Detect as low as 0.5 pg RNA



Wide Dynamic Range – 0.5 – 500 pg RNA



Superior Multiplexing – Plex up to 5 targets per reaction



Eliminate Freezer Storage – up to 9 months stability at room temperature

DESCRIPTION:

Qscript Iyo 1-step is a lyophilized single-reaction reagent optimized for highly sensitive and reproducible one-step RT-qPCR using hydrolysis probes. The reagent contains a hot-start thermo-stable polymerase, a genetically engineered reverse transcriptase as well as other components to ensure higher performance detection of up to 5 targets with maximum sensitivity and specificity. The enhanced stability of the freeze-dried master mix enables convenient shipping and storage at room temperature. The single tube reaction facilitates easy reaction set up while preventing potential cross-contamination.



Convenient setup



Lyophilized
master mix & buffer



Add RNA,
primers/probes
(H₂O optional)

Single tube


No master mix/buffer transfer. Never cross-contaminate.

3.6 Just add RNA, primers and probes.

3.5 Qscript Iyo 1-step lyospheres.




Better sensitivity enables greater detection




3.7 Qscript lyo 1-step demonstrates higher sensitivity and fluorescence signal than a market-leading one-step RT-qPCR product (liquid). The figure shows a LDHA assay with 10-fold serial diluted universal human RNA (50 ng to 0.5 pg) as template. 20 μ l of LDHA primer/probe was used to resuspend each lysosphere, followed by addition of 5 μ l of RNA target. RT-qPCR was carried out on a QuantStudio 7 instrument with the following conditions: 50°C 15 min, 95°C 3 min, 40 cycles: 95°C 15 sec, 60°C 60 sec.

High performance across a range of inputs



3.8 Qscript lyo 1-step was tested with 12 serial dilutions of Human total RNA (4x from 1250 ng to 0.3 pg). 20 μ l of LDHA primer/probe was used to resuspend each lysosphere, followed by addition of 5 μ l of RNA target. RT-qPCR was carried out on a QuantStudio 7 instrument with the following conditions: 50°C 15 min, 95°C 3 min, 40 cycles: 95°C 15 sec, 60°C 60 sec.

Multiplex up to 5 targets per assay



3.9 Qscript lyo 1-step was tested with a 5-plex assay with 7 serial dilutions of Human total RNA (4x from 50 ng to 30 pg). 20 μ l of primers/probes was used to resuspend each lysosphere, followed by addition of 5 μ l of RNA target. RT-qPCR was carried out on a QuantStudio 7 instrument with the following conditions: 50°C 15 min, 95°C 3 min, 40 cycles: 95°C 15 sec, 60°C 60 sec.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
Qscript lyo 1-step - 8 R	95198-008	8 rxns
Qscript lyo 1-step - 24 R	95198-024	24 rxns

Ask about custom sizes and capabilities



UltraPlex One-Step ToughMix

Up to 5-target multiplex, inhibitor-resistant RT-qPCR, maximum yields for superior performance

FEATURES AND BENEFITS:

- 4x concentrated SuperMix reagent supports increased sample input volume, improving flexibility with extremely low yield templates (1 pg total RNA)
- ToughMix reagent technology neutralizes a broad spectrum of PCR inhibitors common in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start
- Temperature stabilized master mix enables convenient setup at ambient room temperature

DESCRIPTION:

UltraPlex One-Step ToughMix is a ready-to-use, single-component SuperMix reagent for One-Step reverse transcription and real-time quantitative PCR (RT-qPCR) of RNA templates using probe-based detection methods. Advanced qScript XLT reverse transcriptase mutant possesses elevated temperature stability and improved template binding affinity increases sensitivity with large, complex RNA targets and delivers highly sensitive

quantification with highly multiplexed RNA detection assays. ToughMix reagent technology ensures robust, reliable performance of highly-multiplexed (>4) RNA detection assays with a wide-range of inhibitory starting materials. This flexible formulation supports miniaturized reaction volumes (droplet PCR) with either fast or standard thermal cycling conditions.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
UltraPlex One-Step ToughMix- 100 R	95166-100	100 x 20 µl rxns
UltraPlex One-Step ToughMix - 500 R	95166-500	500 x 20 µl rxns
UltraPlex One-Step ToughMix - 1000 R	95166-01K	1000 x 20 µl rxns
UltraPlex One-Step ToughMix ROX - 100 R	95167-100	100 x 20 µl rxns
UltraPlex One-Step ToughMix ROX - 500 R	95167-500	500 x 20 µl rxns
UltraPlex One-Step ToughMix Low ROX - 100 R	95168-100	100 x 20 µl rxns
UltraPlex One-Step ToughMix Low ROX - 500 R	95168-500	500 x 20 µl rxns



qScript One-Step SYBR Green RT-qPCR


Sensitive, dye-based RNA quantification with gene specific primers in a single, seamless reaction mixture without opening the tube prior to PCR

FEATURES AND BENEFITS:

- One-step reaction minimizes opportunity for pipetting error
- Robust, specific amplification
- AccuStart hot start mAb technology

DESCRIPTION:

The qScript One-Step SYBR Green RT-qPCR Kit is a convenient and highly sensitive solution for quantitative RT-PCR of RNA templates (RT-qPCR) using SYBR Green I dye detection and gene-specific primers. cDNA synthesis and PCR amplification are carried out in the same tube without opening between procedures. The system has been optimized to deliver maximum RT-PCR efficiency, sensitivity and specificity. The proprietary reaction buffer has been specifically formulated to maximize activities of both reverse transcriptase and Taq DNA polymerase while minimizing the potential for primer-dimer and other non-specific PCR artifacts. The kit is compatible with both fast and standard qPCR cycling protocols. Highly specific amplification is essential for successful RT-qPCR with SYBR Green I technology, since this dye binds to any dsDNA generated during amplification. AccuStart Taq DNA polymerase contains monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step. Upon heat activation at 95°C, the antibodies denature irreversibly, releasing fully active, unmodified Taq DNA polymerase.



3.10 One-Step SYBR Green RT-qPCR with broad dynamic range, high sensitivity and high specificity. A 202 bp fragment of glyceraldehyde-3-phosphate dehydrogenase (GAPD) mRNA was amplified from log-fold serial dilutions of HeLa cell total RNA (100 ng to 0.1 pg). Eight replicate reactions for each RNA quantity, and the no template control (NTC) were carried out in 25 μ l volumes with the qScript One-Step SYBR Green RT-qPCR Kit and 200 nM each GAPD specific primers (PrimerBank ID 7669492a2, Wang, X. and Seed, (2003) NAR 31(24): e154; pp.1-8). Reactions were assembled on ice, transferred to a MyiQ™ real-time detection system (Bio-Rad Laboratories), and incubated for 5 min at 50°C followed by 2 min at 95°C. PCR cycling was for 40 cycles of 3 s, 95°C; 30 s, 60°C. Immediately following PCR cycling the block temperature was ramped from 60°C to 95°C and melt curve data was collected. Panel **A**) Amplification plots and standard curve regression analysis. Panel **B**) Dissociation results (melt curve) for NTC, 0.1 pg and 1 pg reactions.





ORDER INFO

Product Name	Quantabio Catalog Number	Size
qScript One-Step SYBR Green RT-qPCR Kit for iQ - 50 R	95086-050	50 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit for iQ - 200 R	95086-200	200 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit - 50 R	95087-050	50 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit - 200 R	95087-200	200 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit, ROX - 50 R	95088-050	50 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit, ROX - 200 R	95088-200	200 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit, Low ROX - 50 R	95089-050	50 x 50 µl rxns
qScript One-Step SYBR Green RT-qPCR Kit, Low ROX - 200 R	95089-200	200 x 50 µl rxns



qScript XLT One-Step RT-qPCR ToughMix

Robust, inhibitor-resistant RT-qPCR, maximum yields for superior performance


FEATURES AND BENEFITS:

- ToughMix reagent technology neutralizes a broad spectrum of PCR inhibitors common in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start
- Inert AccuVue plate loading dye simplifies reaction setup and provides instant visual cue of reagent addition and mixing
- Supports efficient vortex mixing with proprietary anti-foaming formulation
- Flexible – supports both fast and standard thermal cycling conditions


DESCRIPTION:

qScript XLT One-Step RT-qPCR ToughMix is a ready-to-use, highly sensitive master mix for reverse transcription quantitative PCR (RT-qPCR) of RNA templates using hybridization probe detection chemistries. Advanced qScript XLT reverse transcriptase mutant possesses elevated temperature stability and improved template binding affinity for large complex RNA templates. Supports a wide-range of challenging starting materials with PCR-inhibitor neutralizing ToughMix additives and delivers consistent, reliable assay performance.

Sequential temperature incubations are performed to the same reaction mixture without opening the tube. The proprietary one-step reaction buffer has been specifically formulated to maximize activity of each enzyme while minimizing the potential for primer-dimer and other non-specific PCR artifacts. Inert AccuVue plate loading dye simplifies reaction assembly and provides instant visual confirmation of reagent addition and mixing.



3.11 Broad linear dynamic range, low Limit of Detection.



3.12 Enables performance in the presence of inhibitors.



ORDER INFO

Product Name	Quantabio Catalog Number	Size
qScript XLT One-Step RT-qPCR ToughMix - 100 R	95132-100	100 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix - 500 R	95132-500	500 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix - 2000 R	95132-02K	2000 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix, ROX - 100 R	95133-100	100 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix, ROX - 500 R	95133-500	500 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix, ROX - 2000 R	95133-02K	2000 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix, Low ROX - 100 R	95134-100	100 x 20 µl rxns
qScript XLT One-Step RT-qPCR ToughMix, Low ROX - 500 R	95134-500	500 x 20 µl rxns



qScript One-Step RT-qPCR

Sensitive RNA quantification with probe-based detection chemistries in a single, seamless reaction mixture without opening the tube prior to PCR

FEATURES AND BENEFITS:

- Simplified 2-reagent system supports user-friendly reaction setup at ambient temperature
- Highly sensitive RNA detection with performance engineered, qScript RNase H(+) M-MLV reverse transcriptase mutant
- Superior assay sensitivity and specificity with AccuStart hot start enzyme technology
- Compatible with either fast or standard thermal cycling conditions

DESCRIPTION:

The qScript One-Step RT-qPCR kit is a convenient, 2-component reagent system that supports highly sensitive one-step real-time PCR detection assays of RNA templates (RT-qPCR) and is compatible with all dual-labeled (hydrolysis) probe chemistries. First-strand cDNA synthesis and PCR amplification are carried out in the same tube without opening between procedures. Specifically formulated to maximize activities of both reverse transcriptase and Taq DNA polymerase while minimizing the potential for primer-dimer and other non-specific PCR artifacts; this proprietary one-step formulation has been optimized to deliver maximum RT-PCR efficiency, sensitivity, and specificity,

enabling unbiased co-amplification of low copy transcripts in the presence of higher copy reference genes.

Highly specific amplification is crucial to successful RT-qPCR as non-specific product(s) can compete for amplification of the target sequence and impair PCR efficiency. A key component of this kit is AccuStart Taq DNA polymerase which contains monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step. Upon heat activation at 95°C, the antibodies denature irreversibly, releasing fully active, unmodified Taq DNA polymerase.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
qScript One-Step RT-qPCR Kit - 50 R	95057-050	50 x 50 µl rxns
qScript One-Step RT-qPCR Kit - 200 R	95057-200	200 x 50 µl rxns
qScript One-Step RT-qPCR Kit, ROX - 50 R	95058-050	50 x 50 µl rxns
qScript One-Step RT-qPCR Kit, ROX - 200 R	95058-200	200 x 50 µl rxns
qScript One-Step RT-qPCR Kit, ROX - 1000 R	95058-01K	1000 x 50 µl rxns
qScript One-Step RT-qPCR Kit, Low ROX - 50 R	95059-050	50 x 50 µl rxns
qScript One-Step RT-qPCR Kit, Low ROX - 200 R	95059-200	200 x 50 µl rxns

4.0

PCR & qPCR

4.1

qPCR Instrumentation

Q qPCR Cycler

A faster, smaller, better way to qPCR

FEATURES AND BENEFITS:

-  **Ultra-Fast Data Acquisition** – 35 cycles in 25 minutes*
-  **Unrivaled Performance** – Detect 2-fold expression level differences
-  **Portable & Compact** – 4.5 lbs - transport without ever calibrating
-  **Scalable & Wireless** – Connect up to 10 instruments (48 samples/cycler)
-  **Magnetic Induction** – Eliminate variability vs block-based cyclers

DESCRIPTION:

Q uses a patented magnetic induction technology to rapidly heat samples coupled with fan forced air for cooling to acquire data in as little as 25 minutes. Available in 2 or 4 channel models, the robust optical system acquires all channels simultaneously and allows for running the fastest multiplexed assays.

Q's miniature speaker-size and 4.5 pound weight make it the most portable and versatile qPCR cycler on the market without ever needing to calibrate. Q also provides scalability as each instrument can process up to 48 samples per run and up to 10 Q's can be connected to a single computer wirelessly via bluetooth enabling up to 480 samples to be processed simultaneously.

A key difference is that Q incorporates a unique spinning aluminum rotor providing superior temperature uniformity of $\pm 0.05^{\circ}\text{C}$ versus traditional block-based real time cyclers which rely on multiple peltier heating blocks that can create

edge effects resulting in sample variation. Not only does the data give you superior reproducibility, repeatability but enables detection of 2-fold gene expression level differences as well as identification of difficult class IV SNP's requiring melt temperature resolutions of 0.1°C .



4.1 Q qPCR cycler.



* 25 minute cycle times obtained with fast cycling master mixes and short amplicon assay designs targeting cDNA.



Ultra-Fast Data Acquisition

Generate high quality data, fast!

- Q's speed is the fastest in the industry
- Don't sacrifice on the performance quality of your qPCR
- Completing runs in as little as 25 minutes* is the new standard




4.2 5 point, 2x dilution series of Hepatitis B virus (HBV) cDNA template. Starting amount of 3E+06 copies (n = 4 each); Efficiency = 90% (standard curve method); R² = 0.99; Time to complete run (including melt) = 26 min
*25 minute cycle times obtained with fast cycling master mixes and short amplicon assay designs targeting cDNA

Unrivaled Performance

Detect two-fold differences

- Confidently detect small differences
- High thermal uniformity and reproducibility
- Detect differences within a single cycle



4.3 Manganese superoxide dismutase gene (MnSOD). Eight point, 2x dilution series of human genomic DNA (n = 4 each); Efficiency = 98% (standard curve method); R² = 1.00

ORDER INFO

Product Name

- Q 2-channel qPCR Instrument
- Q 4-channel qPCR Instrument
- Q Tubes & Caps (20 racks/box, total of 960 tubes and caps)

Quantabio Catalog Number

- 95900-2C
- 95900-4C
- 95910-20

Size

- 1 instrument
- 1 instrument
- 1 box

Q Cycler does not require the use of reference dyes.



4.2 Conventional PCR

AccuStart Long Range SuperMix

Superior sensitivity and multiplexing for DNA amplification of long targets

FEATURES AND BENEFITS:


- Amplify +24 kb gDNA and +40 kb lambda DNA
- 4x concentration enables superior sensitivity with low inputs (100 pg)
- Multiplex up to 6 targets with each target up to 6 kb in length
- Stabilized, single-tube SuperMix minimizes pipetting errors and hands-on-time

DESCRIPTION:

The AccuStart Long Range SuperMix is a 4x solution that contains all the components for long range target amplifications, including a blend of two hot-start thermostable DNA polymerases and an optimized buffer. This SuperMix enables routine and easy amplification with high accuracy (>10x Taq) and accommodates targets with broad GC-content (no separate GC buffer needed). This product is also capable of multiplexing and is suitable for End Point PCR, template prep for Sanger Sequencing, NGS, Cloning and HLA typing.


Superior Sensitivity: Improved yields across a range of DNA inputs and target sizes

The AccuStart Long Range SuperMix can amplify DNA inputs as low as 100 pg, across a wide range of target sizes.



LR AccuStart Long Range SuperMix T Takara LR ver2.1
P Promega GoTaq LR

4.4 Comparison of sensitivity and yield. 8.5 kb β -globin fragments were amplified in 50 μ l reaction volumes according to the recommended protocol. Reaction inputs varied from 100 pg – 10 ng. Following a 3 min activation at 95°C; 30 cycles of PCR were performed: 92°C, 30 s; 65°C, 6 min; 72°C, 10 min.




LR AccuStart Long Range SuperMix A NEB LA Master Mix
B KAPA ReadyMix C AccuPrime Taq High Fidelity
D Platinum Taq High Fidelity

4.5 Comparison of yield over fragment length. 8.5 kb β -globin and 24 kb tPA fragments were amplified in 50 μ l reaction volumes according to the recommended protocol. Reaction inputs were 5 ng and 25 ng for the 8 kb and 24 kb fragments, respectively. Following a 3 min activation at 95°C; 27 cycles of PCR were performed: 92°C, 30 s; 65°C, 6 min (8.5 kb), 12 min (24 kb); 72°C, 10 min.



Maximized Multiplexing: Amplify 6 targets up to 6 kb each


The AccuStart Long Range SuperMix can multiplex 6 targets up to 6 kb each, further speeding up experiments, reducing costs and allowing for more data to be derived per run.



4.6 Strong multiplexing capabilities. Multiplex PCR reactions were performed for amplification of 2 kb, 3 kb, 4 kb and 5 kb BRCA 1 targets with 1 kb and 6 kb BRCA 2 targets from 10 ng human genomic DNA template. Reactions were run on MJ Research PTC-200 Thermal Cycler. Following a 2 min activation of 95°C; 27 cycles of PCR were performed: 92°C, 30 s, 65°C, 8 min; 72°C, 10 min. 4 µl of PCR products were analyzed on a 0.7% agarose gel.

Reliable lot-to-lot reproducibility

AccuStart Long Range SuperMix is manufactured in a state of the art facility under an ISO13485 quality system which provides consistent lot-to-lot reproducibility.



- 1 Lot 1
- 2 Lot 2
- 3 Lot 3

4.7 Lot-to-Lot consistency. 8.5 kb β-globin, 17.5 kb and 24 kb tPA fragments were amplified in 50 µl reaction volumes according to the recommended protocol. Reaction inputs were 5 ng, 10 ng and 25 ng for the 8 kb, 17.5 kb and 24 kb fragments, respectively. Following a 3 min activation at 95°C; 27 cycles of PCR were performed: 92°C, 30 s; 65°C, 6 min (8.5 kb), 8 min (17.5 kb), 12 min (24 kb); 72°C, 10 min. 5 µl of PCR products were analyzed on a 0.5% agarose gel with a DNA marker.

ORDER INFO

Product Name

AccuStart Long Range SuperMix - 25 R
AccuStart Long Range SuperMix - 100 R

Quantabio Catalog Number

95199-025
95199-100

Size

25 rxns
100 rxns



repliQa HiFi ToughMix

Superior speed and inhibitor tolerance for DNA amplification requiring high fidelity

FEATURES AND BENEFITS:

- Fidelity of >90x wild type Taq
- 2–3x faster PCR results with extension rates as fast as 1 kb/sec*
- Tough Tested – tolerant to a wide range of PCR inhibitors
- Superior yield and sensitivity
- Amplification of +24 kb gDNA and +40 kb λ DNA


DESCRIPTION:

The repliQa HiFi ToughMix is a 2x, ready-to-use solution that contains all the components for high fidelity PCR amplification, including a genetically modified DNA polymerase coupled with hot start antibodies.

This unique, next generation master mix provides >90x higher fidelity compared to Taq, while reducing time to PCR results by 2–3x. The extreme speed is enabled by extension times as fast as 1–10 kb/sec depending on target length. The enzyme is coupled with the industry leading ToughMix which is tolerant to a wide variety of inhibitors making it suitable for routine PCR, cloning, amplicon sequencing and site directed mutagenesis.

Extreme Speed: 2–3x faster results


repliQa HiFi ToughMix has very fast extension times, ranging from 1–10 kb/sec depending on the fragment size, which can significantly shorten the time to result.



4.8 Comparison of speed. A 2 kb fragment was amplified in 50 μ l reaction volumes according to the recommended protocol. Following a 30 s activation at 98°C; 30 cycles of PCR were performed: 98°C, 10 s; 60°C, 10 s; 68°C, 5–30 s. The thermal cycler had a ramp rate of 5°C/s.

Tough Tested: Tolerant to a wide range of PCR inhibitors


repliQa HiFi ToughMix is able to tolerate a wide range of common PCR inhibitors, allowing for amplification of crude or difficult PCR sample types.



4.9 Strong Inhibitor Resistance. A 2 kb λ DNA template was amplified using each manufacturers recommended cycling conditions with different amounts of inhibitors. The experiment was run in duplicate.

Superior Yield and Sensitivity

repliQa HiFi ToughMix provides higher yield and sensitivity, highlighting the enzyme efficiency. Coupled with extreme amplification speed allows PCR products to be amplified earlier and detected and lower levels.




4.10 Comparison of yield. A gDNA template was amplified with varying GC-content and length targets using each manufacturers recommended cycling conditions. 8 different targets were used. The experiment was run in duplicate.





Superior Yield and Sensitivity

repliQa HiFi ToughMix demonstrates greater efficiency, at almost 90%, enabling higher yields and ultimately better sensitivity.




4.11 Comparison of efficiency. Amplify human genomic DNA template with varying GC-content and length targets using each master mix's recommended cycling conditions. 8 different targets were used. Ran in duplicate. See GC-content slide for more detail.

4.12 Comparison of yield. A gDNA template was amplified with varying GC-content and length targets using each manufacturer's recommended cycling conditions. 8 different targets were used. The experiment was run in duplicate.


Long Amplification

repliQa HiFi ToughMix has the ability to amplify long fragments +24 kb gDNA and +40 kb λ DNA, further proving the versatility of this enzyme.



H: repliQa HiFi ToughMix **P:** NEB Phusion High Fidelity Master Mix **Q:** NEB Q5 High-Fidelity 2x Master Mix **S:** Thermo Invitrogen Platinum SuperFi Master Mix

4.13 Long Range capabilities (gDNA). A range of 3.6 kb, 8.5 kb, 17.5 kb, and 24 kb gDNA templates were amplified with varying GC-content and lengths using each manufacturer's recommended cycling conditions. The experiment was run in duplicate.




4.14 Long Range capabilities (λ DNA). A range of 10 kb, 20 kb, 30 kb, and 40 kb λ DNA templates were amplified with varying GC-content and lengths using each manufacturers recommended cycling conditions. The experiment was run in duplicate.

H: repliQa HiFi ToughMix
 P: NEB Phusion High Fidelity Master Mix
 Q: NEB Q5 High-Fidelity 2x Master Mix
 S: Thermo Invitrogen Platinum SuperFi Master Mix

Consistent GC Tolerance

repliQa HiFi ToughMix is able to amplify varying levels of GC-content targets (32%–70% GC-rich), further enabling superior PCR performance.



4.15 Wide GC-content tolerance range. gDNA templates were amplified with varying GC-content and lengths using each manufacturers recommended cycling conditions. 8 different targets were used. The GC-content varied with 32%/900 base pairs (bp), 34%/800 bp, 38%/900 bp, 40%/1000 bp, 63%/300 bp, 64%/600 bp, 64%/800 bp and 70%/400 bp. The experiment was run in duplicate.

H: repliQa HiFi ToughMix
 P: NEB Phusion High Fidelity Master Mix
 Q: NEB Q5 High-Fidelity 2x Master Mix
 S: Thermo Invitrogen Platinum SuperFi Master Mix

ORDER INFO

Product Name	Quantabio Catalog Number	Size
repliQa HiFi ToughMix - 25	95200-025	25 rxns
repliQa HiFi ToughMix - 100	95200-100	100 rxns
repliQa HiFi ToughMix - 500	95200-500	500 rxns



AccuStart II PCR ToughMix

Robust, reliable PCR assay performance with challenging sample materials or impure templates

FEATURES AND BENEFITS:


- Stabilized 2x PCR SuperMix enables convenient room-temperature setup and is unaffected by repetitive freeze-thaw
- High-yielding, ultrapure modified Taq DNA polymerase delivers robust, reliable duplex assay performance
- Stringent, ultrapure antibody hotstart ensures sensitive and specific target amplification
- Separate electrophoretic mobility dye reduces risk of post-PCR cross contamination with gel electrophoresis

DESCRIPTION:

AccuStart II PCR ToughMix is a 2x concentrated ready-to-use reaction cocktail for robust, general-purpose PCR amplification of DNA templates in the presence of PCR inhibitors. It contains all components, except primers and template. This reagent formulation contains an ultrapure, AccuStart II Taq DNA poly-

merase with stringent antibody hot start to ensure specific and efficient primer extension with convenient reaction assembly at ambient temperature. PCR products generally contain non-templated dA additions and can be cloned using vectors that have a single 3'-overhanging thymine residue on each end.

30 cycle PCR; 1×10^4 copies TcR DNA (1052 bp amplicon)



4.16 Inhibitor Resistance of AccuStart II PCR ToughMix. A 1-kb fragment from 1×10^4 copies of the Tetracyclin resistance gene was amplified in 20 μ l reaction volumes according to the recommended protocol. Reactions were challenged with varying concentrations of different PCR inhibitors as summarized below. Following a 3 min activation at 94°C; PCR was for 30 cycles of: 94°C, 15 s; 60°C, 20 s; 72°C, 1 min. 1/5th of each reaction was analyzed on a 0.1% agarose, 0.5x TBE gel containing 0.25 mg/ml ethidium bromide.

Hemoglobin: 316 ng/ μ l, 100 ng/ μ l, 31.6 ng/ μ l, 10 ng/ μ l, 3.16 ng/ μ l

Humic Acid: 31.6 ng/ μ l, 10 ng/ μ l, 3.16 ng/ μ l, 1 ng/ μ l, 0.316 ng/ μ l


Hematin: 100 μ M, 31.6 μ M, 10 μ M, 3.16 μ M, 1 μ M

Melanin: 10 ng/ μ l, 3.16 ng/ μ l, 1 ng/ μ l, 0.316 ng/ μ l, 0.1 ng/ μ l

C: control reactions without inhibitor; L: 1 Kb Plus DNA Ladder (Invitrogen)



30 cycle PCR; 1 x 10⁴ copies TcR DNA (1052 bp amplicon)



4.17 Inhibitor Resistance of AccuStart II PCR ToughMix: PCR in the presence of polyphenol spike. Varying amounts of a polyphenol-rich plant extract (0.2, 0.06, 0.02, 0.006, or 0.002 µl) were added to 25 µl PCRs containing 10,000 copies of a control template. Amplification was carried out for 30 cycles of: 94°C, 15 s; 60°C, 20 s; 72°C, 1 min. 1/5th of each reaction was analyzed on a 01% agarose, 0.5x TBE gel containing 0.25 mg/ml ethidium bromide. As little as 0.002 µl of the crude plant lysate inhibited control reactions with a conventional PCR master mix (data not shown).

ORDER INFO

Product Name

AccuStart II PCR ToughMix - 100 R
 AccuStart II PCR ToughMix - 800 R
 AccuStart II PCR ToughMix - 4000 R

Quantabio Catalog Number

95142-100
 95142-800
 95142-04K

Size

100 x 25 µl rxns (1 x 1.25 ml)
 800 x 25 µl rxns (8 x 1.25 ml)
 4000 x 25 µl rxns (1 x 50 ml)



AccuStart Taq DNA Polymerase HiFi

Sensitive, precise amplification of large DNAs for high-fidelity downstream applications

FEATURES AND BENEFITS:

- Superior assay sensitivity and specificity with AccuStart enzyme technology – maximum-yielding Taq DNA polymerase mutant with stringent, multi-epitope antibody hot start
- Optimized high-fidelity PCR buffer with proofreading activity supports robust amplification of large PCR products ≤20 kb in length

DESCRIPTION:

AccuStart Taq DNA Polymerase HiFi combines high-yielding mutant Taq DNA polymerase with proof-reading polymerase and ultrapure, monoclonal antibody hot-start activation control. This reagent provides highly-sensitive

and precise target amplification with convenient assembly at ambient temperature. Robust and reliable amplification of large, complex DNA targets up to 20 kb in length.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
AccuStart Taq DNA Polymerase HiFi - 250 U	95085-250	250 units (5 units/μl)
AccuStart Taq DNA Polymerase HiFi - 1000 U	95085-01K	1000 units (5 units/μl)
AccuStart Taq DNA Polymerase HiFi - 5000 U	95085-05K	5000 units (5 units/μl)



AccuStart II GelTrack PCR SuperMix

2x concentrated, ready-to-use reaction cocktail for routine PCR amplification of DNA fragments up to 4 kb

FEATURES AND BENEFITS:

- GelTrack Loading Dye pre-mixed
- High yield, high sensitivity
- Precise amplification – hot-start technology ensures specific and efficient primer extension
- Convenient reaction assembly at room temperature
- Preblended electrophoresis dyes to streamline gel electrophoresis workflows

DESCRIPTION:

AccuStart II GelTrack PCR SuperMix contains all components, except primers and template necessary for robust PCR. It simplifies reaction assembly, improves assay reproducibility, and reduces the risk of contamination. A key component is AccuStart II Taq DNA polymerase which contains monoclonal antibodies that bind to the polymerase and keep it inactive prior to the initial PCR denaturation step. Upon heat activation (1 minute at 94°C), the antibodies denature irreversibly,

releasing fully active, unmodified Taq DNA polymerase. This enables specific and efficient primer extension with the convenience of room temperature reaction assembly.

GelTrack Loading Dye is a mixture of blue and yellow electrophoresis-tracking dyes that migrate at approximately 4kb and 50 bp, and comes pre-mixed with the PCR reagents.

ORDER INFO

Product Name

AccuStart II GelTrack PCR SuperMix - 100 R
AccuStart II GelTrack PCR SuperMix - 500 R
AccuStart II GelTrack PCR SuperMix - 4000 R

Quantabio Catalog Number

95136-100
95136-500
95136-04K

Size

100 x 25 µl rxns (1 x 1.25 ml)
500 x 25 µl rxns (5 x 1.25 ml)
4000 x 25 µl rxns (1 x 50 ml)



AccuStart II PCR SuperMix

Robust, user-friendly 1-tube PCR SuperMix reagents for routine, general purpose PCR

FEATURES AND BENEFITS:

- 1-tube SuperMix reagent minimizes pipetting, simplifies reaction assembly and improves accuracy
- Sensitive, precise DNA amplification with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start

DESCRIPTION:

AccuStart II PCR SuperMix is a 2x concentrated, ready-to-use reaction cocktail for routine PCR amplification of DNA fragments up to 4 kb. 1-tube SuperMix reagent simplifies reaction assembly by minimizing pipetting steps and improving assay reproducibility. Ultrapure AccuStart II Taq DNA polymerase uses

a stringent multi-epitope antibody hot-start that prevents non-specific primer extension prior to heat activation (1 minute at 94°C). The antibodies are irreversibly denatured, releasing a fully active, high-yielding Taq DNA polymerase mutant.

ORDER INFO

Product Name

AccuStart II PCR SuperMix - 100 R
AccuStart II PCR SuperMix - 500 R
AccuStart II PCR SuperMix - 4000 R

Quantabio Catalog Number

95137-100
95137-500
95137-04K

Size

100 x 25 µl rxns (1 x 1.25 ml)
500 x 25 µl rxns (5 x 1.25 ml)
4000 x 25 µl rxns (1 x 50 ml)



AccuStart II Taq DNA Polymerase

High purity, recombinant Taq DNA polymerase preparation with high avidity monoclonal antibodies that bind the polymerase and keep it inactive prior to the initial PCR denaturation step

FEATURES AND BENEFITS:

- Supports specific primer extension with AccuStart technology convenient room temperature reaction assembly

DESCRIPTION:

Upon heat activation (1 minute at 94°C), the antibodies denature irreversibly, releasing fully active, unmodified Taq DNA polymerase. The AccuStart II automatic hot-start enables specific and efficient primer extension in the PCR process with the added convenience of room temperature reaction assembly.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
AccuStart II Taq DNA Polymerase - 250 U	95141-250	250 U (5 U/μl)
AccuStart II Taq DNA Polymerase - 1000 U	95141-01K	1000 U (5 U/μl)
AccuStart II Taq DNA Polymerase - 5000 U	95141-05K	5000 U (5 U/μl)



5PRIME HotMaster Taq DNA Polymerase

Innovative hot-start/cold-stop technology provides highly specific hot start PCR without activation step


FEATURES AND BENEFITS:

- Minimal optimization with self-adjusting magnesium buffer technology
- Novel HotMaster hot-start/cold-stop technology releases full activity at higher temperatures without activation step
- Convenient reaction set-up at ambient temperature

DESCRIPTION:

5PRIME HotMaster Taq DNA Polymerase uses an innovative thermostable Taq inhibitor that releases fully active at high temperatures ensuring precise primer annealing and extension

at optimal temperatures and “cold-stop” after thermal cycling. Ideal for fast thermal cycling protocols.



4.18 Fast PCR Amplification of a 131 bp fragment of the human TNF gene with standard Taq, HotMaster Taq and conventional Hot Start enzymes. PCR protocol: 1 sec 95°C denaturation, 1 sec 55°C annealing, 5 sec 72°C extension. Initial denaturation was 2 min at 95°C prior to PCR or 10 min for the chemically modified enzyme respectively.

ORDER INFO

Product Name

HotMaster Taq DNA Polymerase - 100 U
 HotMaster Taq DNA Polymerase - 1000 U
 HotMaster Taq DNA Polymerase - 5000 U

Quantabio Catalog Number

2200300
 2200320
 2200330

Size

100 U (5 U/μl)
 1000 U (5 U/μl)
 5000 U (5 U/μl)



5PRIME HotMasterMix

Highly specific hot-start PCR with minimal reaction setup

FEATURES AND BENEFITS:

- Optimal results – highly specific amplification
- Minimal handling – ready-to-use mastermix format
- Convenience – storage at 4°C eliminates freeze-thaw cycles

DESCRIPTION:

Highly specific amplification with minimal handling


The 5PRIME HotMasterMix is a ready-to-use reagent mix that offers high reproducibility when processing large numbers of samples. The HotMaster Taq DNA Polymerase, an integral component of the master mix, is designed to reduce or eliminate any non-specific products that may result from mispriming during PCR.

The 5PRIME HotMasterMix (2.5x) is a ready-to-use PCR mix. Only primers and template need to be added to the 2.5x concentrate, thus decreasing the number of time-consuming pipetting steps. This format not only reduces the likelihood of errors and the risk of contamination, but it also increases precision and sample throughput.

The 5PRIME HotMasterMix also contains the self-adjusting Mg^{2+} buffer technology. This formulation adjusts the Mg^{2+} concentration automatically, eliminating the need for optimizing this critical component. The MasterMix does not need to be stored frozen, eliminating the time-consuming thawing process and the resulting reduction in performance.

The 5PRIME HotMasterMix is a 2.5-fold concentrate, resulting in the following final concentrations for a 50 μ l PCR reaction: 1 U Taq DNA Polymerase, 45 mM Cl, 2.5 mM Mg^{2+} , 200 μ M of each dNTP.

Comparison of 5PRIME Taq Enzymes



4.19 Amplification of a 131 bp TNF fragment using different 5PRIME Taq enzymes.

- M 100 bp Marker
- 1 Standard Taq
- 2 HotMaster Taq
- 3 HotMasterMix





ORDER INFO

Product Name	Quantabio Catalog Number	Size
5PRIME HotMasterMix - 100 R	2200400	100 x 50 µl rxns (2 x 1 ml)
5PRIME HotMasterMix - 1000 R	2200410	1000 x 50 µl rxns (20 x 1 ml)

10 mM dNTP Mix

DESCRIPTION:

The 10 mM dNTP Mix is a solution of high purity deoxyribonucleoside triphosphates that has been functionally qualified for real-time quantitative PCR (qPCR). It is suitable for use in conventional end-point PCR, real-time qPCR, first-strand cDNA synthesis, as well as other applications that require dNTP as substrate.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
10 mM dNTP Mix - 200 µl	95062-200	200 µl
10 mM dNTP Mix - 1000 µl	95062-01K	1000 µl



Real-Time Quantitative PCR

4.3

Multiplexed Pre-Amplification

PRODUCT OVERVIEW

	PerfeCra qPCR ToughMix	PerfeCra qPCR ToughMix, UNG	Accustart Genotyping ToughMix	PerfeCra Multiplex qPCR SuperMix
Concentration	2x	2x	2x	2x
Performance	++++	++++	++++	+++
Inhibitor Tolerance	✓	✓	✓	✓
Chemistry	Probe	Probe	Probe	Probe
Sample Type	gDNA, cDNA	gDNA, cDNA	gDNA, cDNA	gDNA, cDNA
Cycling Mode	Standard or Fast	Standard or Fast	Standard or Fast	Standard or Fast
Fast Cycling Compatibility	✓	✓	✓	✓
Multiplex Compatibility	Up to 2 targets	Up to 2 targets	Up to 2 targets	Up to 5 targets
Carryover contamination control	–	Includes heat-labile UNG and a blend of dTTP / dUTP	–	–
Application	Gene Expression	Gene Expression	SNP Genotyping	Gene Expression

	PerfeCra Multiplex qPCR ToughMix	PerfeCra FastMix II	PerfeCra SYBR Green SuperMix	PerfeCra SYBR Green FastMix
Concentration	5x	2x	2x	2x
Performance	++++	+++	++	++
Inhibitor Tolerance	✓	–	–	–
Chemistry	Probe	Probe	SYBR Green	SYBR Green
Sample Type	gDNA, cDNA	gDNA, cDNA	gDNA, cDNA	gDNA, cDNA
Cycling Mode	Standard or Fast	Standard or Fast	Standard only	Standard or Fast
Fast Cycling Compatibility	–	–	–	✓
Multiplex Compatibility	Up to 5 targets	Up to 2 targets	–	–
Carryover contamination control	–	–	–	–
Application	Gene Expression	Gene Expression	Gene Expression, microRNA Expression, ChIP Analysis	Gene Expression Analysis





PerfeCta PreAmp 5x SuperMix

Unbiased pre-amplification of up to 100 DNA targets from as little as 100 pg total cDNA.
Compatible with either probe-based or dye-based qPCR detection chemistries

FEATURES AND BENEFITS:

- Unbiased, linear pre-amplification of up to 100 DNA targets
- 5x concentrated SuperMix maximizes sample input volume with dilute cDNA templates
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot-start
- User-friendly inert AccuVue plate loading dye provides visual confirmation of reagent addition
- Supports efficient vortex mixing

DESCRIPTION:


PerfeCta PreAmp SuperMix is a 5x concentrated, ready-to-use reaction cocktail for unbiased, selected enrichment of target sequences from limiting amounts of starting material for downstream gene expression profiling or targeted re-sequencing. It contains all the necessary components, except primers and templates. The 5x concentrated SuperMix allows addition of higher template volumes when working with low concentration samples, and/or reduced reaction volumes. Inclusion of inert light blue tracer dye helps visualize small reaction volumes and ensure accurate pipetting. PerfeCta PreAmp SuperMix delivers unbiased pre-amplification of up to 100 target

sequences from as little as 100 pg of total cDNA. It is compatible with both TaqMan 5'-nuclease probes or ds-DNA binding dye (i.e. SYBR Green I) qPCR detection chemistries.

A key component of PerfeCta PreAmp SuperMix is an ultra pure, highly processive thermostable DNA polymerase that is combined with high avidity monoclonal antibodies. This proprietary polymerase mix is resistant to PCR inhibitors and provides an extremely stringent automatic hot start, allowing reaction assembly and temporary storage at room temperature prior to pre-amplification.



PreAmp Process Flow



PreAmp Process Flow:

1. Prepare RNA
2. Reverse transcribe RNA
3. Pool assay primers and dilute
4. Perform pre-amplification reaction
5. Dilute PreAmp reaction product
6. Perform individual qPCRs for each pre-amplified gene of interest (GOI)

ORDER INFO

Product Name

PerfeCTa PreAmp 5x SuperMix - 5 R
PerfeCTa PreAmp 5x SuperMix - 40 R

Quantabio Catalog Number

95146-005
95146-040

Size

5 x 50 µl rxns
40 x 50 µl rxns



SYBR Green Detection

PerfeCta SYBR Green SuperMix/FastMix

Sensitive and precise DNA amplification with DNA-intercalating dye based detection chemistry


FEATURES AND BENEFITS:

- 2x concentrated reagents minimize pipetting steps, simplify reaction assembly and improve accuracy
- Superior assay sensitivity and specificity with ultrapure AccuStart enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start
- Supports efficient vortex mixing with proprietary anti-foaming formulation
- FastMix formulation supports both fast and standard thermal cycling conditions
- SuperMix version provides maximum dye concentration for robust optical signal with small amplicons (i.e. microRNA-templated cDNA)


DESCRIPTION:

Specific target amplification is essential for precise target quantification with SYBR Green I technology since this dye binds to all dsDNA generated during amplification. PerfeCta SYBR Green SuperMix and FastMix ensure specific primer extension products with ultra-pure AccuStart hot start technology and proprietary formulation that reduces potential for primer-dimer and other non-specific artifacts.


Single-tube reagents are 2x concentrated ready-to-use reaction cocktails containing all necessary components, except primers and DNA template for quantitative PCR. Proprietary formulation stabilizes SYBR Green I dye to deliver maximum efficiency, sensitivity, and robust fluorescent signal.




4.20 PKCA target amplified from log-fold serial dilutions of qScript synthesized cDNA from HeLa cell total RNA (100 ng – 1 pg) with either PerfeCta SYBR Green SuperMix or Power SYBR Green PCR Master Mix.



4.21 PKCA target amplified from log-fold serial dilutions of qScript synthesized cDNA from HeLa cell total RNA (100 ng – 1 pg) with either PerfeCra SYBR Green SuperMix or SYBR GreenER™ qPCR SuperMix.



4.22 ADAR target amplified from log-fold serial dilutions of qScript synthesized cDNA from HeLa cell total RNA (100 ng – 10 pg) with either PerfeCra SYBR Green FastMix or QuantiFast SYBR Green PCR Kit.



4.23 ADAR target amplified from log-fold serial dilutions of qScript synthesized cDNA from HeLa cell total RNA (100 ng – 10 pg) with either PerfeCra SYBR Green FastMix or iQ™ SYBR Green Supermix.





ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCта SYBR Green FastMix for iQ - 250 R	95071-250	2 x 1.25 ml
PerfeCта SYBR Green FastMix for iQ - 1250 R	95071-012	10 x 1.25 ml
PerfeCта SYBR Green FastMix for iQ - 5000 R	95071-05K	1 x 50 ml
PerfeCта SYBR Green FastMix - 250 R	95072-250	2 x 1.25 ml
PerfeCта SYBR Green FastMix - 1250 R	95072-012	10 x 1.25 ml
PerfeCта SYBR Green FastMix - 5000 R	95072-05K	1 x 50 ml
PerfeCта SYBR Green FastMix, ROX - 250 R	95073-250	2 x 1.25 ml
PerfeCта SYBR Green FastMix, ROX - 1250 R	95073-012	10 x 1.25 ml
PerfeCта SYBR Green FastMix, ROX - 5000 R	95073-05K	1 x 50 ml
PerfeCта SYBR Green FastMix, Low ROX - 250 R	95074-250	2 x 1.25 ml
PerfeCта SYBR Green FastMix, Low ROX - 1250 R	95074-012	10 x 1.25 ml
PerfeCта SYBR Green FastMix, Low ROX - 5000 R	95074-05K	1 x 50 ml
PerfeCта SYBR Green SuperMix for iQ - 250 R	95053-100	2 x 1.25 ml
PerfeCта SYBR Green SuperMix for iQ - 1250 R	95053-500	10 x 1.25 ml
PerfeCта SYBR Green SuperMix for iQ - 5000 R	95053-02K	1 x 50 ml
PerfeCта SYBR Green SuperMix - 250 R	95054-100	2 x 1.25 ml
PerfeCта SYBR Green SuperMix - 1250 R	95054-500	10 x 1.25 ml
PerfeCта SYBR Green SuperMix - 5000 R	95054-02K	1 x 50 ml
PerfeCта SYBR Green SuperMix, ROX - 250 R	95055-100	2 x 1.25 ml
PerfeCта SYBR Green SuperMix, ROX - 1250 R	95055-500	10 x 1.25 ml
PerfeCта SYBR Green SuperMix, ROX - 5000 R	95055-02K	1 x 50 ml
PerfeCта SYBR Green SuperMix, Low ROX - 250 R	95056-100	2 x 1.25 ml
PerfeCта SYBR Green SuperMix, Low ROX - 1250 R	95056-500	10 x 1.25 ml
PerfeCта SYBR Green SuperMix, Low ROX - 5000 R	95056-02K	1 x 50 ml



Probe-based Detection

PerfeCta qPCR ToughMix

Robust, inhibitor-resistant probed-based qPCR

FEATURES AND BENEFITS:

- 2x concentrated reagents minimize pipetting steps, simplify reaction assembly and improve accuracy
- ToughMix reagent technology neutralizes a broad spectrum of PCR inhibitors common in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start
- Easy-to-use 2x concentrated master mixes with AccuVue plate loading dye and optimized passive reference dye for simplified reaction setup
- Supports efficient vortex mixing with proprietary anti-foaming technology

DESCRIPTION:

PerfeCta qPCR ToughMix is a 2x concentrated qPCR SuperMix ready-to-use reaction cocktail for PCR amplification of DNA templates that overcomes a broad spectrum of PCR inhibitors often encountered in environmental specimens, plant tissues or animal tissues. This proprietary polymerase mix provides maximum sensitivity and PCR efficiency with all dual-label (hydrolysis) probe-based detection chemistries and stringent hot-start activation control allowing reaction assembly and pre-run storage at ambient room temperature prior to thermal cycling. Inert AccuVue plate loading dye is compatible with either white or clear PCR plates and helps to minimize pipette error and provides visual confirmation of reagent addition. UNG containing versions are blended with Uracil N-glycosylase to eliminate potential post-PCR carryover contamination associated with routine testing workflows.


Inhibitor	Common sources	Reagent performance	
		Competitor	PerfeCta ToughMix
Polyphenols	Plant extracts	–	✓
Humic acids	Soil Plant tissues	–	✓
Hematin	Dried bloods Blood spots	–	✓
Hemoglobin	Blood	✓	✓
Polysaccharides	Feces Plant tissues	–	✓
Melanin	Hair Skin	–	✓






Effect of PCR Inhibitors on qPCR of MYC cDNA

PerfeCta qPCR ToughMix




Slope	-3.365	-3.267	-3.420
R ²	-0.9989	-0.9970	-0.9988
Int.	27.57	27.42	28.08
PCR Eff.	98.2%	102.4%	96.1%

Path-ID™ qPCR Master Mix



Slope	-3.376	-3.935	-3.258
R ²	-0.9968	-0.9953	-0.9989
Int.	27.07	33.93	29.32
PCR Eff.	97.7%	79.5%	102.7%

TaqMan Environmental Master Mix 2.0



Slope	-3.409	-3.616	-3.480
R ²	-0.9968	-0.9965	-0.9989
Int.	28.14	30.77	28.92
PCR Eff.	96.5%	89.0%	93.8%

4.24 Serial dilutions of qScript cDNA. cDNA alone, cDNA + 1 µg hemoglobin, cDNA + 10 ng/µl humic acid (100 ng/rxn). 10 µl reactions; Roche LC480; 384-well optimal cycling for TaqMan reagents: 95°C, 10 min; followed by 45 cycles of 95°C, 15 s; 60°C, 60 s. 0.5x MYC (FAM-MGB) TaqMan Gene Expression Assay.



ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCта qPCR ToughMix - 250 R	95112-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix - 1250 R	95112-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix - 5000 R	95112-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCта qPCR ToughMix, ROX - 250 R	95113-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix, ROX - 1250 R	95113-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix, ROX - 5000 R	95113-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCта qPCR ToughMix, Low ROX - 250 R	95114-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix, Low ROX - 1250 R	95114-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix, Low ROX - 5000 R	95114-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCта qPCR ToughMix, UNG - 250 R	95138-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG - 1250 R	95138-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG - 5000 R	95138-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCта qPCR ToughMix, UNG, ROX - 250 R	95139-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG, ROX - 1250 R	95139-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG, ROX - 5000 R	95139-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCта qPCR ToughMix, UNG, Low ROX - 250 R	95140-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG, Low ROX - 1250 R	95140-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCта qPCR ToughMix, UNG, Low ROX - 5000 R	95140-05K	5000 x 20 µl rxns (1 x 50 ml)



PerfeCta FastMix II

Convenient, 1-tube reagent solution supports robust and reliable probe-based DNA detection with fast or standard thermal cycling conditions

FEATURES AND BENEFITS:

- 2x concentrated reagents minimize pipetting steps, simplify reaction assembly and improve accuracy
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot-start
- Inert AccuVue plate loading dye simplifies reaction setup and provides instant visual cue of reagent addition and mixing
- Supports efficient vortex mixing with proprietary anti-foaming technology
- Flexible – optimized for both fast and standard thermal cycling conditions

DESCRIPTION:

PerfeCta FastMix II is an advanced 2x concentrated qPCR SuperMix reagent for both fast and conventional PCR cycling protocols or instruments. It is a versatile and robust solution that provides the superior sensitivity and high PCR efficiency and is compatible with all dual-label (hydrolysis) probe chemistries. The kit is provided as a 2x concentrated ready-to-use reaction cocktail that contains all required reaction components, except primers, probe(s), and q DNA template. Inert AccuVue

plate loading dye helps to minimize pipette error and provides visual confirmation of thorough mixing.

A key component of the kit is an ultrapure, highly-processive thermostable DNA polymerase that is free of detectable *E. coli* DNA. PerfeCta FastMix II is ideal for demanding qPCR applications such as bacterial pathogen detection where residual host DNA in typical recombinant enzyme preparations can limit assay sensitivity and obscure detection of low copy samples.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCta qPCR FastMix II - 250 R	95118-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCta qPCR FastMix II - 1250 R	95118-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCta qPCR FastMix II - 5000 R	95118-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCta qPCR FastMix II, ROX - 250 R	95119-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCta qPCR FastMix II, ROX - 1250 R	95119-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCta qPCR FastMix II, ROX - 5000 R	95119-05K	5000 x 20 µl rxns (1 x 50 ml)
PerfeCta qPCR FastMix II, Low ROX - 250 R	95120-250	250 x 20 µl rxns (2 x 1.25 ml)
PerfeCta qPCR FastMix II, Low ROX - 1250 R	95120-012	1250 x 20 µl rxns (10 x 1.25 ml)
PerfeCta qPCR FastMix II, Low ROX - 5000 R	95120-05K	5000 x 20 µl rxns (1 x 50 ml)



PerfeCta MultiPlex qPCR ToughMix

Advanced 1-tube SuperMix optimized to support highly multiplexed DNA amplifications in miniaturized reaction volumes and withstand a broad spectrum of PCR inhibitors


FEATURES AND BENEFITS:

- 1-tube SuperMix minimizes pipetting steps, simplifying reaction assembly and improving accuracy
- 5x concentrated reagent provides greater sensitivity and more flexibility with dilute DNA samples
- ToughMix reagent technology neutralizes a broad spectrum of PCR inhibitors common in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hotstart
- Supports efficient vortex mixing with proprietary anti-foaming additives

DESCRIPTION:

Suppression of low copy amplicons by high copy reference targets during multiplex co-amplification skews the apparent representation and quantification of low copy target sequences. PerfeCta MultiPlex qPCR ToughMix transcends these limitations by enabling sensitive, broad linear dynamic detection range with co-amplification of four abundant (10^6) targets. PerfeCta MultiPlex qPCR ToughMix enables multiplex qPCR assay performance comparable to single-plex qPCR probe assays without the need to rigorously titrate primer concentration.

PerfeCta MultiPlex qPCR ToughMix is a 5x concentrated, ready-to-use reaction cocktail for real-time quantitative PCR (qPCR) with ToughMix reagent technology that neutralizes a broad spectrum of PCR inhibitors that compromise assay performance with crude extracts, clinical specimens, plants, soil, and environmental or complex food matrices. A key component of the kit is an ultrapure, highly-processive thermostable DNA polymerase that is free of detectable *E. coli* DNA. PerfeCta MultiPlex qPCR ToughMix is ideal for demanding qPCR applications such as bacterial pathogen detection where residual host DNA in typical recombinant enzyme preparations can limit assay sensitivity and obscure detection of low copy samples.



4.25 High efficiency, high sensitivity multiplex qPCR results with PerfeCta MultiPlex qPCR ToughMix. Log-fold serial dilutions (10 to 1 E7 copies) of a plasmid containing the GAPDH gene, as well as no template controls, were amplified with PerfeCta MultiPlex qPCR ToughMix as either a single-plex qPCR, or a 4-target multiplexed qPCR that contained 1 E8 copies of 3 additional plasmid DNAs (ACTB, IL1beta, and TUBA). Quadruplicate reactions for each input quantity were carried out in 25 μ l volumes with 300 nM each primer and 150 nM each probe. Dual-labeled probes with non-fluorescent quenchers were from Biosearch Technologies. GAPDH was detected using a FAM-BHQ1 probe. ACTB; CAL Fluor Orange 560-BHQ1; IL1beta: CAL Fluor Red 610-BHQ2; TUBA: Quasar 670 – BHQ3. Single-plex qPCRs only contained the GAPDH primers and probe. Cycling was performed on a Bio-Rad CFX with the following protocol. 95°C, 2 min; followed by 40 cycles of 95°C, 10 s; 58°C, 90 s. RFU data were exported to Excel, averaged for each replicate reaction series, and plotted.





ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCта MultiPlex qPCR ToughMix - 250 R	95147-250	250 x 25 µl rxns (1 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix - 1000 R	95147-01K	1000 x 25 µl rxns (4 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix - 5000 R	95147-05K	5000 x 25 µl rxns (1 x 25 ml)
PerfeCта MultiPlex qPCR ToughMix, ROX - 250 R	95148-250	250 x 25 µl rxns (1 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix, ROX - 1000 R	95148-01K	1000 x 25 µl rxns (4 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix, ROX - 5000 R	95148-05K	5000 x 25 µl rxns (1 x 25 ml)
PerfeCта MultiPlex qPCR ToughMix, Low ROX - 250 R	95149-250	250 x 25 µl rxns (1 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix, Low ROX - 1000 R	95149-01K	1000 x 25 µl rxns (4 x 1.25 ml)
PerfeCта MultiPlex qPCR ToughMix, Low ROX - 5000 R	95149-05K	5000 x 25 µl rxns (1 x 25 ml)



PerfeCta MultiPlex qPCR SuperMix

Sensitive and robust multiplex qPCR assay performance with user-friendly 1-tube reagents containing pre-blended passive reference dye

FEATURES AND BENEFITS:

- Robust assay performance for highly-multiplexed DNA quantification assays
- 2x concentrated reagents minimize pipetting steps, simplify reaction assembly and improve accuracy
- Superior assay sensitivity and specificity with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hot start
- Broad linear detection range with highly-multiplexed qPCR assays
- Supports efficient vortex mixing with proprietary anti-foaming additives

DESCRIPTION:

Suppression of low abundance targets by high abundance reference targets during co-amplification is a common problem in multiplex PCR in which individual assay sensitivity can be significantly compromised. PerfeCta Multiplex qPCR SuperMix delivers assay performance with exceptionally broad, linear detection and limit-of-detection (LOD) sensitivity.

Multiplexed qPCR performance comparable to single-plex assay performance without the need for rigorous titration of individual primer assays is achievable with this reagent.

PerfeCta Multiplex qPCR SuperMix is a 2x concentrated, ready-to-use reaction cocktail that contains all the necessary

components except: primers, probe(s), and DNA template for highly-multiplexed, real-time quantitative PCR. This reagent formulation pushes the boundary of multiplex qPCR by enabling unbiased amplification of up to 5 targets in a single amplification reaction. A key component of the kit is an ultrapure, highly-processive thermostable DNA polymerase that is free of detectable *E. coli* DNA. PerfeCta MultiPlex qPCR SuperMix is ideal for demanding qPCR applications such as bacterial pathogen detection where residual host DNA in typical recombinant enzyme preparations can limit assay sensitivity and obscure detection of low copy samples.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCta MultiPlex qPCR SuperMix - 50 R*	95063-050	50 x 50 µl rxns (1 x 1.25 ml)
PerfeCta MultiPlex qPCR SuperMix - 200 R*	95063-200	200 x 50 µl rxns (4 x 1.25 ml)
PerfeCta MultiPlex qPCR SuperMix - 1000 R*	95063-01K	1000 x 50 µl rxns (1 x 25 ml)
PerfeCta MultiPlex qPCR SuperMix, Low ROX - 50 R	95108-050	50 x 50 µl rxns (1 x 1.25 ml)
PerfeCta MultiPlex qPCR SuperMix, Low ROX - 200 R	95108-200	200 x 50 µl rxns (4 x 1.25 ml)
PerfeCta MultiPlex qPCR SuperMix, Low ROX - 1000 R	95108-01K	1000 x 50 µl rxns (1 x 25 ml)

* Contains separate tube of 50x ROX and 50x Low Rox



4.4

Cloning

repliQa HiFi ToughMix

Superior speed and inhibitor tolerance for DNA amplification requiring high fidelity

FEATURES AND BENEFITS:


- Fidelity of >90x wild type Taq
- 2–3x faster PCR results with extension rates as fast as 1 kb/sec*
- Tough Tested – tolerant to a wide range of PCR inhibitors
- Superior yield and sensitivity
- Amplification of +24 kb gDNA and +40 kb λ DNA

DESCRIPTION:

The repliQa HiFi ToughMix is a 2x, ready-to-use solution that contains all the components for high fidelity PCR amplification, including a genetically modified DNA polymerase coupled with hot start antibodies. This unique, next generation master mix provides >90x higher fidelity compared to Taq, while reducing time to PCR results by 2–3x. The extreme speed is enabled by extension times as fast as 1–10 kb/sec depending on target length. The enzyme is coupled with the industry leading ToughMix which is tolerant to a wide variety of inhibitors making it suitable for routine PCR, cloning, amplicon sequencing and site directed mutagenesis.

Extreme Speed: 2–3x faster results


repliQa HiFi ToughMix has very fast extension times, ranging from 1–10 kb/sec depending on the fragment size, which can significantly shorten the time to result.



4.8 Comparison of speed. A 2 kb fragment was amplified in 50 μ l reaction volumes according to the recommended protocol. Following a 30 s activation at 98°C; 30 cycles of PCR were performed: 98°C, 10 s; 60°C, 10 s; 68°C, 5–30 s. The thermal cycler had a ramp rate of 5°C/s.

Tough Tested: Tolerant to a wide range of PCR inhibitors


repliQa HiFi ToughMix is able to tolerate a wide range of common PCR inhibitors, allowing for amplification of crude or difficult PCR sample types.



4.27 Strong Inhibitor Resistance. A 2 kb λ DNA template was amplified using each manufacturers recommended cycling conditions with different amounts of inhibitors. The experiment was run in duplicate.

Superior Yield and Sensitivity

repliQa HiFi ToughMix provides higher yield and sensitivity, highlighting the enzyme efficiency. Coupled with extreme amplification speed allows PCR products to be amplified earlier and detected and lower levels.




4.10 Comparison of yield. A gDNA template was amplified with varying GC-content and length targets using each manufacturers recommended cycling conditions. 8 different targets were used. The experiment was run in duplicate.





Superior Yield and Sensitivity

repliQa HiFi ToughMix demonstrates greater efficiency, at almost 90%, enabling higher yields and ultimately better sensitivity.




4.11 Comparison of efficiency. Amplify human genomic DNA template with varying GC-content and length targets using each master mix's recommended cycling conditions. 8 different targets were used. Ran in duplicate. See GC-content slide for more detail.

4.12 Comparison of yield. A gDNA template was amplified with varying GC-content and length targets using each manufacturer's recommended cycling conditions. 8 different targets were used. The experiment was run in duplicate.


Long Amplification

repliQa HiFi ToughMix has the ability to amplify long fragments +24 kb gDNA and +40 kb λ DNA, further proving the versatility of this enzyme.



H: repliQa HiFi ToughMix **P:** NEB Phusion High Fidelity Master Mix **Q:** NEB Q5 High-Fidelity 2x Master Mix **S:** Thermo Invitrogen Platinum SuperFi Master Mix

4.13 Long Range capabilities (gDNA). A range of 3.6 kb, 8.5 kb, 17.5 kb, and 24 kb gDNA templates were amplified with varying GC-content and lengths using each manufacturer's recommended cycling conditions. The experiment was run in duplicate.




4.14 Long Range capabilities (λ DNA). A range of 10 kb, 20 kb, 30 kb, and 40 kb λ DNA templates were amplified with varying GC-content and lengths using each manufacturers recommended cycling conditions. The experiment was run in duplicate.

H: repliQa HiFi ToughMix
 P: NEB Phusion High Fidelity Master Mix
 Q: NEB Q5 High-Fidelity 2x Master Mix
 S: Thermo Invitrogen Platinum SuperFi Master Mix

Consistent GC Tolerance

repliQa HiFi ToughMix is able to amplify varying levels of GC-content targets (32%–70% GC-rich), further enabling superior PCR performance.



4.15 Wide GC-content tolerance range. gDNA templates were amplified with varying GC-content and lengths using each manufacturers recommended cycling conditions. 8 different targets were used. The GC-content varied with 32%/900 base pairs (bp), 34%/800 bp, 38%/900 bp, 40%/1000 bp, 63%/300 bp, 64%/600 bp, 64%/800 bp and 70%/400 bp. The experiment was run in duplicate.

H: repliQa HiFi ToughMix
 P: NEB Phusion High Fidelity Master Mix
 Q: NEB Q5 High-Fidelity 2x Master Mix
 S: Thermo Invitrogen Platinum SuperFi Master Mix

ORDER INFO

Product Name	Quantabio Catalog Number	Size
repliQa HiFi ToughMix - 25	95200-025	25 rxns
repliQa HiFi ToughMix - 100	95200-100	100 rxns
repliQa HiFi ToughMix - 500	95200-500	500 rxns



repliQa HiFi Assembly Mix

Seamless assembly of multiple DNA fragments for high efficiency cloning

FEATURES & BENEFITS:

- Formulated to increase number of transformants
- Assemble up to 6 fragment inserts without the need for restriction enzymes
- Flexibility in design with 10x master mix enabling assembly of low concentration DNA samples
- Eliminates dilution step resulting in easier workflow and 1 hr assembly
- Includes DpnI to reduce background when using plasmid templates for PCR


DESCRIPTION:

The repliQa HiFi Assembly Mix simplifies the construction of recombinant DNA through seamless assembly of multiple DNA fragments in a single, isothermal reaction.

Similar in principle to the Gibson Assembly® Method, the high efficiency repliQa HiFi Assembly Mix is ideal for a range of genetic engineering applications including:

- Routine molecular cloning
- Site-directed mutagenesis
- Synthetic biology
- Construction of libraries for directed evolution studies

The concentrated (10x), two component format allows flexibility in design of assembly reactions and compatibility with less concentrated DNA samples. The repliQa Mix has been optimized for use with a total input quantity of DNA fragments in the range of 0.03 to 0.5 pmols. Assembly of up to six DNA fragments is recommended, though the repliQa Mix has been used successfully for more complex assemblies. The mix supports assembly of multiple DNA fragments in a single 1-hr reaction.






Speed up your workflow and increase transformation efficiency

repliQa HiFi Assembly Mix increases transformation efficiency without the need for diluting or purifying the assembly reaction prior to transformation of competent cells, resulting in less hands-on time and faster workflows.

Low DNA amounts can be used and will efficiently generate high amounts of transformants.

Dilution of the assembly reaction not required – high transformation efficiency, less hands-on time




4.34 Three DNA fragments containing 23 bp overlaps were generated by PCR, DpnI treated and purified. The three fragments, 4.2 kb, 3.1 kb, and 400 bp in size, were combined in a 1:1.4:5 molar ratio. Total DNA quantities used are indicated (x-axis) and reacted at 50 °C for 60 minutes according to the protocol. One microliter of the undiluted assembled products or one microliter of a 4-fold dilution of the assembled products was used to transform 30 µl of chemically competent cells.

Assemble larger number of fragments

repliQa HiFi Assembly Mix allows for large inserts to be cloned at a very high efficiency resulting in significantly more positive clones eliminating the need to repeat experiments due to erroneous or insufficient clones.

High transformation efficiency results in greater number of fragments



4.35 PCR fragments containing 30 bp overlaps were DpnI treated, purified, and assembled according to the protocol. Reactions contained the indicated number of DNA fragments (0.1 pmol each) and were incubated at 50 °C for 60 minutes. 1 µl of the reactions were used to transform 30 µl of chemically competent cells.

ORDER INFO

Product Name

repliQa HiFi Assembly Mix, 10 rxns
repliQa HiFi Assembly Mix, 50 rxns

Quantabio Catalog Number

95190-010
95190-050

Size

10 rxns
50 rxns

AccuStart Genotyping ToughMix

AccuStart Genotyping ToughMix enables probe-based genetic analysis (SNP detection and allelic discrimination) directly from crude extracts, DBS punches, plant tissue and clinical specimens

FEATURES AND BENEFITS:


- Optimized buffer chemistry destabilizes single base-pair mismatch probes, providing superior allelic discrimination and improved cluster separation for critical, single-nucleotide polymorphism (SNP) detection assays
- Sensitive, precise detection with ultrapure AccuStart II enzyme technology – maximum-yielding Taq DNA polymerase mutant controlled by stringent, multi-epitope antibody hotstart
- ToughMix reagent technology neutralizes a broad spectrum of PCR inhibitors common in plant and animal tissues, environmental sources, clinical specimens and complex food matrices
- Easy-to-use 2x concentrated SuperMix with AccuVue plate loading dye and pre-blended passive reference dye simplifies reaction setup
- Supports efficient vortex mixing with proprietary anti-foaming technology

DESCRIPTION:

Genotyping ToughMix is a 1-tube qPCR SuperMix reagent compatible with all dual-label (hydrolysis) probe chemistries for both fast and conventional PCR cycling protocols or instruments. This proprietary formulation has been rigorously optimized to destabilize single base-pair mismatches to ensure precise allelic discrimination and cluster separation with SNP detection


assays. The reagent is provided as a 2x concentrated ready-to-use reaction cocktail that contains all required reaction components, except primers, probe(s), and DNA template. Inert AccuVue plate loading dye helps to minimize pipette error and provides visual confirmation of thorough mixing.

ToughMix vs Competitor



5.1 Comparison to conventional master mixes AccuStart Genotyping ToughMix stands up to the challenge where other genotyping master mixes fall apart. ToughMix can be used with clean templates where it generates higher fluorescent signal and tighter clusters than the leading competitors.

Influence of PCR inhibitor



5.2 Comparison to conventional master mixes. In the presence of a common PCR inhibitor, humic acid (50 ng/μl), the competitors system is completely shut down while ToughMix delivers robust, accurate results



ORDER INFO

Product Name	Quantabio Catalog Number	Size
Genotyping ToughMix - 250 R	95115-250	250 x 20 μ l rxns (2 x 1.25 ml)
Genotyping ToughMix - 1250 R	95115-012	1250 x 20 μ l rxns (10 x 1.25 ml)
Genotyping ToughMix - 5000 R	95115-05K	5000 x 20 μ l rxns (1 x 50 ml)
Genotyping ToughMix, ROX - 250 R	95116-250	250 x 20 μ l rxns (2 x 1.25 ml)
Genotyping ToughMix, ROX - 1250 R	95116-012	1250 x 20 μ l rxns (10 x 1.25 ml)
Genotyping ToughMix, ROX - 5000 R	95116-05K	5000 x 20 μ l rxns (1 x 50 ml)
Genotyping ToughMix, Low ROX - 250 R	95117-250	250 x 20 μ l rxns (2 x 1.25 ml)
Genotyping ToughMix, Low ROX - 1250 R	95117-012	1250 x 20 μ l rxns (10 x 1.25 ml)
Genotyping ToughMix, Low ROX - 5000 R	95117-05K	5000 x 20 μ l rxns (1 x 50 ml)



AccuStart II PCR Genotyping Kit

Completely reagent-based system enables reliable PCR genotyping with minimal pipetting skill

FEATURES AND BENEFITS:

- Premixed electrophoresis mobility loading dye reduces chances for post-PCR cross contamination
- Stabilized 2x PCR SuperMix enables convenient room-temperature setup
- High-yielding, ultrapure modified Taq DNA polymerase delivers robust, reliable duplex assay performance
- Stringent, ultrapure antibody hotstart ensures sensitive and specific target amplification
- Flexible protocol delivers rapid results in as little as 10 minutes


DESCRIPTION:

The AccuStart II Genotyping Kit is a complete reagent kit designed to support conventional, end-point PCR-based screening of transgenic animal models commonly used in life science research and is validated for use with mouse, fish, or insect tissue specimens. It combines a rapid, 2-component DNA extraction reagent with a user-friendly 2x concentrated PCR SuperMix with loading dye for seamless gel electrophoresis analysis. qPCR-grade genomic DNA template is obtained with

minimal extraction volumes ($\leq 100 \mu\text{l}$) and can be carried out in ≤ 30 -minutes on a standard PCR thermal cycler.

Contents

- Extracta DNA Prep for PCR - Tissue, 2 x 25 ml Extraction Reagent and 2 x 25 ml Stabilization Buffer
- AccuStart II GelTrack PCR SuperMix (2x)
- 5 x 1.25 ml of 2x reaction mixture containing optimized concentrations of MgCl_2 , dNTPs, AccuStart II Taq DNA Polymerase, AccuStart Taq antibodies, reaction buffer, stabilizers and gel loading dyes. Individual components can be reordered separately.



5.3 Two mouse tail snips (2.5 mm) were extracted according to the recommended conditions for each kit. The volume of each extract was brought to 300 μl and diluted 1/20 with TE buffer. 5 μl of diluted extract was used in a 25 μl PCR reactions.

ORDER INFO

Product Name

AccuStart II Genotyping Kit - 100 R
AccuStart II Genotyping Kit - 500 R

Quantabio Catalog Number

95135-100
95135-500

Size

100 x 25 μl rxns
500 x 25 μl rxns



AccuMelt HRM SuperMix

AccuMelt HRM SuperMix maximizes differences in melt temperature and curve shape to allow discrimination of DNA sequence differences amongst different samples

FEATURES AND BENEFITS:

- See sequence differences clearly – robust amplification ensures sufficient yield of products to generate discrete melt curves
- Accurate genotype calling – comparable or better performance than TaqMan Genotyping
- Work with rare or precious samples – large range of template inputs possible
- Specificity – works with lower Mg^{2+} concentration than other systems thus enhancing assay accuracy


DESCRIPTION:

AccuMelt HRM SuperMix



AccuMelt HRM SuperMix is a ready-to-use 2x concentrated hot-start PCR mix containing SYTO 9™ green fluorescent DNA-binding dye.

Superior Resolution of Genotypes

SNP Genotyping is a useful application for HRM and illustrates the capabilities of AccuMelt HRM SuperMix. Genotypes are readily identified based on unique melting profiles depending on a sample's sequence (Figure 1). Furthermore, AccuMelt HRM SuperMix gives superior resolution of difficult genotypes when compared to the leading competitor's mix based on greater T_m differences observed for A → T transversions (Figure 2).



5.4 High resolution melting analysis of a model SNP system with a single A,C,G, or T variant base. AccuMelt HRM SuperMix readily resolves each genotype and T_m differences are easily visualized in normalized melting curve plots (Roche, Lightcycler 480).



5.5 Effect of T,A,C, or G variant base on T_m in a model HRM SNP system with either AccuMelt HRM SuperMix (Panel A), or a competitor's SYTO 9 dye master mix (Panel B). Plots of averaged melt peaks normalized to maximum signal for each system.






Comparison to TaqMan Genotyping

TaqMan Genotyping has been used successfully in SNP analysis and other allelic discrimination applications. This widely adopted standard in genotyping was used as a benchmark to assess the utility of HRM with our SuperMix. AccuMelt HRM

was determined to be just as effective as TaqMan Genotyping in SNP analysis and was even able to call the genotype for a difficult sample which the TaqMan assay could not resolve.

Comparison to TaqMan Genotyping




5.6 Accuracy of HRM genotyping with AccuMelt HRM SuperMix was evaluated by comparison to TaqMan detection of the G>A rs1801133 SNP in the MTHFR gene (Panel A). HRM normalized melting curves (Panel B). TaqMan allelic discrimination plots. TaqMan failed to resolve Sample D3 (labeled as “UNKN”) which was typed as a heterozygote by HRM.

Robust Amplification

Consistent robust amplification is critical to accuracy in HRM analysis. AccuMelt HRM SuperMix will drive all PCR amplifications to plateau regardless of the quantity of template input. This ensures accurate results regardless of the quantity of DNA available.

5.7 High yield, high efficiency PCR with AccuMelt HRM SuperMix. Real-time PCR of GAPDH was amplified from log-fold serial dilutions of qScript synthesized cDNA from HeLa cell total RNA (10 ng to 0.1 pg) was carried out with either a leading SYBR Green Master Mix or AccuMelt HRM SuperMix using the following cycling conditions: 95°C, 20 s; followed by 45 cycles of: 95°C, 3 s; 60°C, 20 s. Averaged plots for quadruplicate reactions for each input quantity are shown.



ORDER INFO

Product Name

AccuMelt HRM SuperMix - 250 R
AccuMelt HRM SuperMix - 1250 R

Quantabio Catalog Number

95103-250
95103-012

Size

250 x 20 µl rxns (2 x 1.25 ml)
1250 x 20 µl rxns (10 x 1.25 ml)

Sensitive, precise quantification of microRNAs using optimized reagent solutions and economical dye-based qPCR detection

The qScript microRNA cDNA Synthesis Kit is optimized to provide highly-sensitive reverse transcription of small, non-coding RNA and is compatible with either total RNA samples or miRNA-enriched template. This optimized kit provides reagent components for polyadenylation of small non-coding RNAs and first-strand cDNA synthesis with qScript reverse transcriptase. A novel sequence is incorporated by the oligo(dT) primer that can be utilized for qPCR quantification using the PerfeCta Universal microRNA Primer and miR-specific forward assay primer (PerfeCta miR Assay). The kit comes complete with a SNORD44 Human Positive Control Primer and extra Poly(A) Tailing Buffer and MicroRNA cDNA Reaction Mix to enable minus poly(A) polymerase control reactions to enable precise control for signal contributions from precursor pri-microRNA.

PerfeCta SYBR Green SuperMixes contain ultrapure, high-yielding Taq with stringent antibody hotstart and maximum SYBR Green dye content to enable sensitive and precise DNA amplification of small PCR amplicons (e.g. microRNA templated cDNA).

qScript microRNA cDNA Synthesis Kit

Reverse transcribe small, single-stranded RNAs with superior sensitivity

FEATURES AND BENEFITS:

- Log-fold greater sensitivity than inefficient stem-loop priming methods.
- Broad linear dynamic detection range across a range of RNA inputs (10 pg – 1 µg total RNA).
- Incorporate novel adapter sequence that can be utilized for qPCR quantification using PerfeCta Universal microRNA Primer

DESCRIPTION:


The qScript microRNA cDNA Synthesis Kit is an optimized reagent system that reverse transcribes small single-stranded RNA into 5'-labeled cDNA using either total RNA or miRNA enriched samples.

Single-stranded RNA is first polyadenylated by poly(A) polymerase before reverse transcription into universal cDNA using high performance qScript RT with a proprietary adapter oligo(dT) primer. The universal cDNA template enables simple and cost-effective microRNA profiling when used together with wet-lab validated PerfeCta microRNA Assays, PerfeCta Universal PCR Primer and PerfeCta SYBR Green SuperMix.


This complete kit includes positive (human) control primer, SNORD44, that can be used to quantitate ubiquitously expressed small nucleolar RNA. In addition, the kit contains 20% extra poly(A) polymerase reaction buffer and microRNA cDNA Reaction Mix to perform (-) poly(A) polymerase and (-) reverse transcriptase control reactions.



Quantabio's qScript microRNA system




miR-1 in heart




6.1 Detection of Rare miRNAs. microRNAs Quantabio's microRNA profiling system provides linear detection and quantification of miRNAs across total-RNA input levels spanning six orders of magnitude. This means miRNAs will be detected even when tissue is scarce or the miRNAs are rare. The Quantabio miRNA system will detect low copy miRNAs more reliably than other systems and due to the minus-PAP control you will have confidence in the results.




miR-124a



miR-27a



Let 7a



6.2 Superior Sensitivity. Quantabio's miRNA Profiling System yields superior results when compared to other leading miRNA quantification systems. The qScript microRNA cDNA Synthesis Kit, PerfeCta microRNA-specific assays and PerfeCta SYBR Green SuperMix form an integrated system that yields industry-leading results.

ORDER INFO

Product Name

qScript microRNA cDNA Synthesis Kit - 25 R
 qScript microRNA cDNA Synthesis Kit - 100 R

Quantabio Catalog Number

95107-025
 95107-100

Size

25 x 20 µl rxns
 100 x 20 µl rxns



Primer design using the qScript miRNA Quantification System

Basic Steps for Primer Design

1. Convert miRNA sequence to a DNA sequence
2. Append the reverse complement of the oligo-dT adapter primer (sequence provided below) to its 3' end.
3. Using primer design software or web-based tool, design the miRNA-specific FORWARD PRIMER that is compatible and T_m -balanced with the REVERSE PRIMER (universal primer, sequence provided below).

Example – human miR-193P

1.

miRNA sequence:

5'-UGUGCAAUUCUAUGCAAACUGA-3'



Convert to DNA sequence

5'-TGTGCAAATCTATGCAAACCTGA-3'

2.

Oligo dT Adapter Primer:

5'-GCATAGACCTGAATGGCGGTAAGGGTGTGGTAGCGGAGACATTTTTTTTTTTTTTTTTTTT-3'



Make reverse complement

5'-AAAAAAAAAAAAAAAAAATGTCTCGCTACCACACCCTTACCGCCATTAGGTCTATGC-3'

Template created by appending the RC of adapter primer to miRNA (DNA seq):

5'-TGTGCAAATCTATGCAAACCTGAAAAAAAAAAAAAAAAAAATGTCTCGCTACCACACCCTTACCGCCATTAGGTCTATGC-3'

3.

Design FORWARD PRIMER

TGTGCAAATCTATGCAAACCTGA →

5'-TGTGCAAATCTATGCAAACCTGAAAAAAAAAAAAAAAAAAATGTCTCGCTACCACACCCTTACCGCCATTAGGTCTATGC-3'

← ATGGCGGTAAGTCCAGATACG

Select last 21 bases as REVERSE PRIMER (Universal Primer)

When designing the forward primer, specific to your miRNA, restrict the search to 1-25 bases. In the example above the software (Oligo7) determined the optimal primer sequence (green) to be in effect the miRNA sequence. In other cases, depending on T_m balancing and compatibility with the universal primer, it may be longer or shorter.

Once you have designed the your miR-specific forward primer, you can order it through the oligonucleotide vendor of your choice. The reverse (universal) primer is included with the qScript microRNA cDNA synthesis kit. The experimental protocol for using the primers is available at www.quantabio.com/resources, along with protocols for the qScript microRNA cDNA synthesis kit and PerfeCta SYBR Green SuperMix.

Developed Assays

A list of primers for more than 1500 miRNAs have already been designed and are available for download at www.quantabio.com/products/microna-profiling. If you are a customer who previously purchased our PerfeCta microRNA Assays, you will find the sequence for the assays of interest in this document. Use the ctrl-F function within excel to search by ID, accession number, or miRNA sequence.

sparQ DNA Frag & Library Prep Kit

Integrated enzymatic fragmentation and library prep with unrivaled speed and performance

FEATURES AND BENEFITS:

- High quality libraries in 2.5 hours from 1 ng – 1 µg of input DNA
- Tunable and reproducible fragmentation size range
- Simple, convenient 2-step workflow with minimal hands-on time
- Novel chemistry and high-fidelity amplification minimizing bias
- Superior sequence coverage uniformity and low duplication rate

DESCRIPTION:


The sparQ DNA Frag & Library Prep Kit is optimized for enzymatic fragmentation of DNA and streamlined construction of high quality libraries for sequencing on Illumina® NGS platforms. The simple, convenient 2-step workflow can be completed in 2.5 hours with minimal hands-on time and accommodates DNA input amounts from 1 ng to 1 µg.

Quantabio's engineered DNA frag and polishing enzymes work in concert to generate fragment sizes that are tunable and reproducible based on reaction time. The DNA fragmentation

and polishing reactions are combined in a single step to generate 5'-phosphorylated and 3'-dA-tailed fragments. This minimizes over fragmentation while greatly simplifies the library prep workflow. Subsequent high efficiency ligation of adapters is performed in the same tube without an intervening cleanup step. If library amplification is required, the HiFi PCR Master Mix and Primer Mix ensure even amplification with minimal bias.

This kit is compatible with single-indexed, or dual-indexed Y-shaped adapters routinely used in library construction.

Streamlined workflow (1 ng – 1 µg input DNA)



7.1 The streamlined workflow utilizes a proprietary enzyme mix that combines fragmentation and DNA polishing in a single step to simplify library construction. High efficiency adapter ligation and cleanup are performed in the same tube, followed by an optional amplification step using HiFi PCR Master Mix and Primer Mix.




Tunable & reproducible fragmentation


The sparQ DNA Frag & Library Prep kit is designed to produce consistent and reproducible fragments that are tunable to application-specific sizes. The fragmentation profile closely resembles Covaris mechanical shearing. Flexible input DNA amounts ranging from 1 ng – 1 µg can be accommodated. The

single-tube enzyme mix fragments DNA and then automatically proceeds to the DNA polishing reaction, thus minimizing potential over fragmentation. Guidelines of incubation time and expected size based on input amount are provided below.

Fragmentation Time Course



Fragmentation Tuning Guide



7.2 sparQ DNA Frag & Library Prep Kit is tunable to the desired fragment size. 100 ng Human gDNA was subjected to fragmentation with a series of incubation time points (3 – 18 min). After fragmentation, DNA samples were purified and then visualized using an Agilent High Sensitivity DNA Kit.


7.3 Guidelines for tuning fragmentation size. If input DNA falls between values displayed on the graph, an estimate can be used for optimizing fragmentation times.

Superior library prep efficiency and yields

The novel and optimized chemistry used in the sparQ DNA Frag & Library Prep Kit coupled with proprietary enzyme mix lead to better sensitivity and higher library yields. PCR-free workflows are enabled for 100 ng of input DNA. For applications requiring amplification, the HiFi PCR Master Mix and Primer Mix

allow researchers to achieve target concentration with very few cycles thereby reducing PCR-derived artifacts. Ultimately, precious samples can be saved for additional applications when necessary.

Workflow Yield Comparison



7.4 sparQ DNA Frag & Library Prep Kit shows significantly higher NGS library preparation efficiency. Libraries with 300 bp average DNA fragments from 100 ng of gDNA Coriell NA12878 were prepared using sparQ DNA Frag & Library Prep Kit and a commercial kit from Supplier K. Manufacturers' manuals were carefully followed. Amplified libraries (5 PCR cycles) were quantified by Qubit fluorometric quantitation method.




Uniform coverage across a wide GC-spectrum


The sparQ chemistry enables high quality library construction with even coverage across a board GC-spectrum including challenging GC- and AT-rich regions. Reproducible and uniform genome coverage is achieved independent of input DNA amounts, comparable to coverage obtained by mechanical

shearing workflows. The sparQ DNA Frag & Library Prep Kit ensures similar total coverage depth for the majority of genomic targets, thus reducing the need for additional sequencing, resulting in less sequencing per sample and lower total sequencing costs.


Genome Coverage Analysis (1 ng input DNA)



Coverage Distribution Analysis (100 ng input DNA)



Genome Coverage Analysis (100 ng input DNA)



7.5 Library prepared using sparQ DNA Frag & Library Prep Kit resulted in uniform coverage across the wide range of GC-content. Libraries were prepared using different DNA fragmentation and library preparation kits with 1 ng or 100 ng of microbial gDNA followed by sequencing on Illumina MiSeq. 2 million reads from each tested library were down-sampled and analyzed. Coverage uniformity against GC-content resulting from different DNA fragmentation and library preparation kits were compared by plotting normalized coverage for a wide GC-content. Libraries prepared using PCR-free workflow of sparQ DNA Frag & Library Prep Kit with 100 ng of microbial genomic DNA shows similar high performance as a typical amplified library prepared by mechanical shearing method.





High quality sequencing metrics with low duplication rates

Excellent sequencing metrics – high mapping percentage and low duplication artifacts – are achieved with the sparQ DNA Frag & Library Prep Kit, ensuring the greatest return on sequencing investments.

	Fragmentation	1 ng input DNA		100 ng input DNA	
		Mapped reads	Duplication	Mapped reads	Duplication
sparQ	Enzymatic	91.9%	0.07%	94.5%	0.04%
Supplier K	Enzymatic	92.4%	0.08%	93.5%	0.03%
Supplier I	Tagmentation	93.8%	0.28%	–	–
Supplier C	Mechanical	93.0%	0.09%	93.6%	0.03%

sparQ DNA Frag & Library Prep Kit generates high quality DNA libraries with minimal duplication artifacts. Libraries were prepared with 1 ng and 100 ng of microbial genomic DNA, amplified for 12 and 6 cycles respectively, and subsequently sequenced on Illumina MiSeq. Each library was down-sampled to 2 million reads (150 bp paired-end reads) and aligned to a reference genome with only unique alignments reported.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
sparQ DNA Frag & Library Prep Kit - 24	95194-024	24 rxns
sparQ DNA Frag & Library Prep Kit - 96	95194-096	96 rxns
sparQ DNA Library Prep Kit - 24	95191-024	24 rxns
sparQ DNA Library Prep Kit - 96	95191-096	96 rxns
Related Products		
sparQ HiFi PCR Master Mix	95192-050	50 rxns (1 x 1.25 ml)
sparQ HiFi PCR Master Mix	95192-250	250 rxns (5 x 1.25 ml)
sparQ PureMag Beads - 5 ml	95196-005	5 ml
sparQ PureMag Beads - 60 ml	95196-060	60 ml
sparQ PureMag Beads - 450 ml	95196-450	450 ml
sparQ Fast Library Quant Kit (for Q) - 50 R	95197-050	50 rxns
sparQ Fast Library Quant Kit (for Q) - 500 R	95197-500	500 rxns
PerfeCta® NGS Quantification Kit - Illumina - 500 R	95154-500	500 x 20 µl rxns
PerfeCta NGS Quantification Kit - Illumina, ROX - 500 R	95155-500	500 x 20 µl rxns
PerfeCta NGS Quantification Kit - Illumina, Low ROX - 500 R	95156-500	500 x 20 µl rxns



sparQ DNA Library Prep Kit

Streamlined, versatile single-tube solution for high quality library prep

FEATURES AND BENEFITS:


- Fast, easy single-tube solution completes library prep in 2.5 hours
- Suitable for a wide range of input amounts from as low as 250 pg
- Optimized chemistry ensuring superior library prep sensitivity and efficiency
- Higher library yields compared to other library prep kits
- High efficiency enables PCR-free workflow from 100 ng input

DESCRIPTION:

The sparQ DNA Library Prep Kit is optimized for the rapid construction of DNA libraries from fragmented double-stranded DNA for sequencing on Illumina® NGS platforms. The simplified protocol speeds up library prep to 2.5 hours with minimal hands-on time and accommodates DNA input amounts from 250 pg to 1 µg. DNA polishing reactions are streamlined into a single step to convert fragmented DNA into 5'-phosphorylated and

3'-dA-tailed DNA fragments. This is followed by high efficiency adapter ligation in the same tube. PCR-free workflows are enabled from 100 ng of starting material. If library amplification is required, the HiFi PCR Master Mix and Primer Mix ensure even amplification with minimal bias. The kit is compatible with input amounts from 250 pg to 1 µg DNA and multiple sample types.

Streamlined workflow (250 pg – 1 µg input DNA)



7.6 The streamlined workflow combines DNA polishing and adapter ligation in a single tube for rapid construction of libraries from fragmented DNA. An optional step using the HiFi PCR Master Mix and Primer Mix ensures even library amplification with minimal bias.







Higher library yields and consistent library prep efficiency

The quality of a library prep depends heavily on the efficiency and sensitivity of the enzymes involved in DNA polishing and adapter ligation steps. The sparQ enzymes have been engineered for optimal sensitivity and efficiency, supporting the

construction of adapter-ligated libraries from a broad range of input DNA from as little as 250 pg. The unique proprietary cocktail of enzymes ensures exceptional library yields – 50% more than other commercial library prep kits.



7.7 sparQ DNA Library Prep Kit produces high quality libraries from a broad range of DNA inputs with significantly higher yields. Libraries were prepared with Covaris-sheared human genomic DNA (250 bp average size) using kit manufacturers' instructions. Amplified libraries (6 PCR cycles for 100 ng input DNA and 13 PCR cycles for 1 ng input DNA) were quantified with Qubit fluorometric method.




7.8 sparQ DNA Library prep Kit resulted in consistent library prep efficiency across a broad range of sample inputs. Libraries were prepared from Covaris-sheared human genomic DNA with sparQ DNA Library Prep Kit without library amplification. Preamplified libraries were quantified with qPCR-based method.

Even coverage across a broad range of GC-content

The sparQ DNA Library Prep Kit enables the construction of high quality libraries with uniform coverage across a wide range of GC-content. For applications requiring amplification, the HiFi

PCR Master mix is formulated to increase library yields while reducing the number of cycles required to create a sequencing-ready library, thereby minimizing PCR-derived artifacts.




7.9 Amplified libraries were prepared from 100 ng of microbial genomic DNA and subsequently sequenced on Illumina MiSeq. 2 million reads from each tested library were down-sampled and analyzed. Coverage uniformity for different library preparation kits were compared by plotting normalized coverage for both extreme AT-rich regions (8%–20% GC-content) and GC-rich regions (75%–88% GC-content).




Coverage of amplified library prepared with sparQ DNA Library Prep Kit closely resembles coverage of PCR-free workflows, both indicating good representation of GC- and AT-rich regions in the final library. Even coverage ensures greater sequencing depth or multiplexing capabilities.


F. nucleatum (27% GC) + *E. coli* (50% GC)
+ *B. pertussis* (67% GC)



F. nucleatum (27% GC)



B. pertussis (67% GC)



■ Amplified library
 ■ PCR-free library
 ■ GC distribution

7.10 Library amplification with sparQ HiFi PCR Master Mix contained in the sparQ DNA Library Prep Kit resulted in uniform coverage across the wide range of GC-content. Libraries were prepared by using sparQ DNA Library Prep Kit with 100 ng input DNA. Coverage depth against GC-content of libraries amplified by sparQ HiFi PCR Master Mix (orange) were compared to corresponding libraries without amplification (dark blue: PCR-free library). GC-content distribution of targeted genomes is indicated by the gray line.

Improved sequencing results

Whether using an amplified or a PCR-free workflow, sparQ DNA Library Prep Kit produces industry leading sequencing results as determined by the high number of reads mapping back to the reference genome and minimal duplication rates.

	Library amplification	1 ng input DNA		100 ng input DNA	
		Mapped reads	Duplication	Mapped reads	Duplication
sparQ	With amplification	94.3%	0.07%	95.5%	0.04%
Supplier K		95.0%	0.09%	95.6%	0.04%
Supplier N		94.9%	0.07%	95.4%	0.03%
sparQ	PCR-free			95.6%	0.03%
Supplier K				95.3%	0.02%
Supplier N				95.1%	0.02%

sparQ DNA Library Prep Kit generates high quality DNA libraries with minimal duplication rates. Libraries were prepared with 1 ng and 100 ng of microbial genomic DNA and subsequently sequenced on Illumina MiSeq. Each library was down-sampled to 2 million reads (150 bp paired-end reads) and aligned to a reference genome with only unique alignments reported.





ORDER INFO

Product Name	Quantabio Catalog Number	Size
sparQ DNA Library Prep Kit - 24	95191-024	24 rxns
sparQ DNA Library Prep Kit - 96	95191-096	96 rxns

Related Products

sparQ HiFi PCR Master Mix	95192-050	50 rxns (1 x 1.25 ml)
sparQ HiFi PCR Master Mix	95192-250	250 rxns (5 x 1.25 ml)
sparQ PureMag Beads - 5 ml	95196-005	5 ml
sparQ PureMag Beads - 60 ml	95196-060	60 ml
sparQ PureMag Beads - 450 ml	95196-450	450 ml
sparQ Fast Library Quant Kit (for Q) - 50 R	95197-050	50 rxns
sparQ Fast Library Quant Kit (for Q) - 500 R	95197-500	500 rxns
PerfeCta® NGS Quantification Kit - Illumina - 500 R	95154-500	500 x 20 µl rxns
PerfeCta NGS Quantification Kit - Illumina, ROX - 500 R	95155-500	500 x 20 µl rxns
PerfeCta NGS Quantification Kit - Illumina, Low ROX - 500 R	95156-500	500 x 20 µl rxns



sparQ HiFi PCR Master Mix

High-fidelity library amplification while maintaining even coverage

FEATURES AND BENEFITS:

- HiFi DNA polymerase engineered to minimize amplification bias
- Increased amplification efficiency resulting in higher yields
- Uniform coverage across challenging AT- and GC-rich regions
- Robust amplification from input DNA as low as 250 pg

DESCRIPTION:

The sparQ HiFi PCR Master Mix is a high-fidelity, high-efficiency PCR master mix for NGS workflows requiring DNA library amplification prior to sequencing. The included primer mix allows amplification of DNA libraries flanked by adapters containing the P5 and P7 sequences required for Illumina® sequencing platforms. The hot-start, proofreading DNA polymerase used


in the sparQ HiFi PCR Master Mix is specifically engineered to improve library amplification efficiency while reducing PCR-derived artifacts, resulting in higher library yields and better coverage uniformity. This kit supports low DNA input from 250 pg and efficient amplification of AT- and GC-rich regions with minimal bias.

Higher library amplification efficiency

Specially designed for sensitive, high efficiency library amplification from a broad range of DNA input, the sparQ HiFi PCR Master Mix minimizes the number of amplification cycles


needed to achieve the threshold required for sequencing. The result is >45% higher library yields while reducing PCR-derived artifacts.

Library Yield Analysis



7.11 Library amplification with sparQ HiFi PCR Master Mix resulted in higher yields. Libraries were prepared from Covaris-sheared human genomic DNA with sparQ DNA library prep kit prior to library amplification. Pre-amplified libraries were then amplified using sparQ HiFi PCR Master Mix (orange) or equivalent kit from Supplier K (blue) and Supplier N (red) with identical PCR cycle numbers (6 cycles for 100 ng input DNA, 14 cycles for 1 ng input DNA, and 16 cycles for 250 pg input DNA). Amplified libraries were quantified with Qubit fluorometric method and qPCR-based quantification method (data not shown).

DNA Libraries from 250 pg Input DNA




7.12 sparQ HiFi PCR Master Mix demonstrates high efficiency library amplification from low input. The fragment size distribution and the quality of the amplified DNA libraries from 250 pg input DNA were analyzed using a high sensitivity DNA analysis kit on the Agilent BioAnalyzer. Libraries were amplified using sparQ HiFi PCR Master Mix (orange) or equivalent kit from Supplier K (blue) and Supplier N (red) with identical amplification cycle numbers (16 cycles for 250 pg input DNA).




Superior coverage uniformity

Libraries amplified by sparQ HiFi PCR Master Mix provide uniform coverage across a broad range of GC-content, similar to corresponding libraries without PCR. Even coverage ensures greater sequencing depth or multiplexing capabilities.


F. nucleatum (27% GC) + *E. coli* (50% GC)
+ *B. pertussis* (67% GC)



F. nucleatum (27% GC)



B. pertussis (67% GC)



7.13 Consistent coverage over a broad range of GC-content with sparQ HiFi PCR Master Mix. Libraries were prepared by using sparQ DNA Library Prep Kit with 100 ng input DNA. Coverage depth against GC-content of libraries amplified by sparQ HiFi PCR Master Mix (red) were compared to corresponding libraries without amplification (dark blue: PCR-free library). GC-content distribution of targeted genomes is indicated by the gray line.

— Amplified library — PCR-free library — GC distribution

ORDER INFO

Product Name

sparQ HiFi PCR Master Mix - 50 R
sparQ HiFi PCR Master Mix - 250 R

Quantabio Catalog Number

95192-050
95192-250

Size

50 rxns (1 x 1.25 ml)
250 rxns (5 x 1.25 ml)

Related Products

sparQ DNA Frag & Library Prep Kit - 24 R
sparQ DNA Frag & Library Prep Kit - 96 R
sparQ DNA Library Prep Kit - 24 R
sparQ DNA Library Prep Kit - 96 R
sparQ PureMag Beads - 5 ml
sparQ PureMag Beads - 60 ml
sparQ PureMag Beads - 450 ml
sparQ Fast Library Quant Kit (for Q) - 50 R
sparQ Fast Library Quant Kit (for Q) - 500 R

95194-024 24 rxns
95194-096 96 rxns
95191-024 24 rxns
95191-096 96 rxns
95196-005 5 ml
95196-060 60 ml
95196-450 450 ml
95197-050 50 rxns
95197-500 500 rxns



sparQ PureMag Beads

Fast, reliable DNA purification & size selection for NGS workflows

FEATURES AND BENEFITS:


- High recovery of DNA fragments greater than 100 bp
- Efficient removal of unwanted components from adapter ligation and PCR reactions
- Consistent single or double-sided size selection
- Seamless integration into existing NGS workflows with little or no protocol change

DESCRIPTION:

sparQ PureMag Beads is a fast and reliable nucleic acid purification system for reaction cleanup and size selection in Next Generation Sequencing (NGS) workflows. Based on the reversible nucleic acid-binding properties of magnetic beads, this product can be used to quickly remove primers, primer-dimers, unincorporated nucleotides, salts, adapters and adapter-dimers from NGS library prep reactions to improve downstream sequencing performance.

sparQ PureMag Beads allows excellent recovery of fragments greater than 100 bp without centrifugation or filtration. Consistent and reliable size selection can be achieved by simply adjusting the beads to sample ratio. This product is designed for both manual and automated processing, allowing seamless integration into existing workflows.


sparQ PureMag Beads workflow



7.14 Double-sided size selection is used to remove smaller and larger fragments from either side of the desired region. The fragment size can be easily adjusted to suit the application by manipulating the sparQ PureMag Beads to DNA volumetric ratio.





Bioanalyzer trace of fragmented human genomic DNA pre- and post double-sided size selection



sparQ PureMag Beads DNA ratio	Targeted size range (bp)	Peak average (bp)
0.7x, 0.9x	200–400	301
0.6x, 0.8x	250–500	377
0.5x, 0.7x	300–700	464


7.15 Electropherogram of fragmented human genomic DNA pre- and post double-sided size selection. Different sparQ PureMag Beads to DNA ratios were used to achieve various targeted size range.

Highly reproducible purification across a range of inputs



7.16 Highly reproducible DNA library profiles were achieved using different lots of sparQ PureMag Beads and a broad range of input amount. Libraries were prepared with sparQ DNA Library Prep Kit from 100 ng and 1 ng of fragmented microbial genomic DNA. sparQ PureMag Beads were used post adapter ligation and PCR amplification to effectively remove adapter-dimers and primer-dimers.

Efficient recovery of DNA equivalent to AMPure XP



7.17 sparQ PureMag Beads show equivalent performance to AMPure XP for DNA purification. 50 bp DNA ladder was purified with sparQ PureMag Beads and AMPure XP at different beads to DNA ratios and analyzed on 2% agarose gel.



ORDER INFO

Product Name	Quantabio Catalog Number	Size
sparQ PureMag Beads - 5 ml	95196-005	5 ml
sparQ PureMag Beads - 60 ml	95196-060	60 ml
sparQ PureMag Beads - 450 ml	95196-450	450 ml



sparQ Fast Library Quant Kit (for Q)

Fastest qPCR-based library quantification in 40 minutes

FEATURES AND BENEFITS:

- Faster time to results – 50% shorter run time than traditional cycling protocols
- Accurate, reliable quantification of NGS libraries of various sizes and GC-content
- High amplification efficiency across a wide linear dynamic range
- Stabilized, ready-to-use sparQ Fast Mastermix to reduce pipetting steps
- Superior run to run uniformity ensuring highly precise measurements



40 min



DESCRIPTION:

sparQ Fast Library Quant Kit provides rapid and accurate quantification of libraries prepared for sequencing on Illumina® NGS platforms. Accurate quantification of the number of amplifiable library molecules prior to loading onto a flow cell is a critical step in the NGS workflow and it ensures optimal cluster genera-

tion and cost-effective use of sequencing capacity. The sparQ Fast Library Quant Kit uses real-time quantitative PCR (qPCR) to specifically quantify the number of library molecules that possess the appropriate adapter tag at each end.

Accurate library quantification in 40 minutes

This kit is optimized for the Q qPCR instrument which uses a magnetic induction technology to rapidly heat samples coupled with fan forced air for cooling to acquire data more rapidly.

The combination of the sparQ Fast Library Quant Kit and the Q instrument enables fast cycling, reducing qPCR run time by 50% compared to traditional cycling protocols.



7.18 Comparison of average qPCR run time for library quantification. sparQ Fast Library Quant Kit uses fast cycling protocol, allowing results to be achieved in 40 minutes versus 1 hour and 20 minutes with the traditional NGS Library Quant Kit.


Complete library quantification solution with unmatched convenience

sparQ Fast Library Quant Kit contains six stabilized, pre-diluted DNA standards, ready-to-use 1.25x mastermix pre-mixed with primer sets containing Illumina P5 and P7 sequences, and an

optimized buffer for diluting NGS library samples. This unique formulation minimizes pipetting steps and ensures precise qPCR results.



sparQ Fast Library Quant Kit workflow



7.19 Illustration of sparQ Fast Library Quant Kit workflow. Reactions are prepared by simply adding standard or diluted library sample. Optimized protocols with fast cycling condition are provided for both 10 μ l or 20 μ l reaction volumes.

qPCR as the most accurate method for library quantification

Real-time quantitative PCR is the most sensitive and precise method for quantifying adapter-ligated DNA molecules. Other methods based on spectrophotometry, fluorometry, or microfluidic electrophoresis (e.g. nanodrop, Qubit, or Bioanalyzer) are acceptable for estimating the appropriate dilutions to use for library quantification. These methods, however, are prone to variabilities and inaccuracies due to factors such as sensitivity to


contaminants or measurement of unsequenceable fragments. The qPCR-based sparQ Fast Library Quant Kit measures only library DNA fragments containing the appropriate adapter sequences on both ends and thus serves as the most accurate method for library quantification, which in turn facilitates optimal loading onto sequencing flow cells.



7.20 Comparison of commonly used methods for measuring library concentration.




High amplification efficiency across a wide linear dynamic range



7.21 sparQ Fast Library Quant Kit provides high amplification efficiency across a wide linear dynamic range. A 10-fold dilution series was prepared from libraries of low (27%) and high (67%) GC-content and amplified under fast conditions on the Q using the sparQ Fast Mastermix. The slopes of the Cq vs Log (concentration) plots and the individual sample reactions measured by the LinRegPCR algorithm indicated superb amplification efficiencies.


Outstanding repeatability of multiple runs



7.22 Multiple runs using sparQ Fast Library Quant Kit resulted in similar quantification values. Two NGS library preparations with high and low GC-contents were amplified and quantified in five distinct runs. Plots of the Normalized fluorescence vs Cycles and average quantification values show the high repeatability of measurements with the sparQ Fast Library Quant Kit under fast cycling on the Q.




Equivalent performance with 50% faster run time



7.23 Results from sparQ and Roche KAPA Library Quant Kits were highly correlated. Concentrations of six different diluted libraries were determined using either the sparQ Fast Library Quant Kit on Q or the Roche KAPA Library Quantification Kit on Bio-Rad CFX following the manufacturer's recommended protocol. Run times, including melt curves, were 40 minutes for sparQ and 80 minutes for KAPA.

Lot-to-lot consistency of sparQ DNA Standards



7.24 sparQ Fast Library Quant Kits are manufactured with high lot-to-lot consistency. Concentrations of diluted libraries with low GC (library 1), high GC (library 2), or balanced GC-content (library 3) were determined using 3 different lots of sparQ DNA Standards. Each library sample was tested in quadruplicate reactions with each lot of sparQ DNA Standards. Standard deviations of average quantification values were all <0.13 pM.

ORDER INFO

Product Name

sparQ Fast Library Quant Kit - 50
sparQ Fast Library Quant Kit - 500

Quantabio Catalog Number

95197-050
95197-500

Size*

50 rxns
500 rxns

* Based on 20 μ l reaction volume.



PerfeCtra NGS Quantification Kit

Real-Time qPCR is the most sensitive and precise method for quantifying adapter-ligated DNA molecules

PerfeCtra NGS Quantification kits provide a comprehensive reagent solution in a single easy-to-use kit

FEATURES AND BENEFITS:

- Precise quantification of adapter-tagged library molecules
- Accurate and sensitive method for NGS library quantification
- Stabilized, prediluted standards for convenient use
- Consistency across a broad range of samples

DESCRIPTION:

Accurate quantification of the number of amplifiable library molecules is the most critical step in the NGS workflow in obtaining high quality read data with next-generation sequencing technologies. The PerfeCtra NGS Quantification Kit uses real-time PCR to specifically quantify library molecules that possess the appropriate adapter tag at each end. These are the suitable template molecules for Bridge PCR used for Illumina NGS platforms. PerfeCtra NGS Quantification Kits simplify the library quantification process by providing stabilized, pre-diluted standards, pre-qualified primer sets, and an optimized dilution buffer for your NGS library samples. This minimizes pipetting errors and ensures reproducible and precise qPCR results, even with dilute samples. The robust qPCR performance of PerfeCtra SYBR® Green SuperMix provides accurate quantification of NGS libraries with varying fragment sizes or GC content.

A common problem with some NGS library quantification protocols is the use of DNA standards that are too concentrated and generate qPCR data that are outside of the linear dynamic range for many qPCR instruments. Improper baseline settings result in compressions between the highest concentrated DNA standards that in turn give rise to inflated PCR efficiencies and inaccurate library quantification results.

The NGS DNA standards supplied with the PerfeCtra NGS Quantification Kits have been carefully selected to avoid these artifacts and produce NGS library standard curves with exceptionally high linear regression correlation coefficients.

Illumina


The DNA standard for Illumina NGS platforms generates a 426 bp amplicon (48.8% GC). Primer sequences correspond to the "P5" and "P7" primer sequences for Illumina sequencing libraries:

Illumina forward primer: 5'-AAT GAT ACG GCG ACC ACC GA-3'


Illumina reverse primer: 5'-CAA GCA GAA GAC GGC ATA CGA-3'



A



B



7.25 PerfeCta NGS Quantification Kit performance data. qPCR amplification of each of the five supplied DNA standards for Illumina NGS libraries (**panel A**) were carried out with the supplied primer sets (300 nM final concentration) and PerfeCta SYBR Green SuperMix in 20 μ l reaction volumes on a Bio-Rad CFX-96. Reactions were incubated for 5 min at 95°C followed by 35 cycles of: 95°C, 10 s; 60°C, 20 s; 45 s, 72°C. Real-time fluorescence data was collected and analyzed at completion of the 72°C extension step. After completion of PCR, a dissociation (melt) curve was performed to verify amplification of a single specific product (**panels B**).

ORDER INFO

Product Name	Quantabio Catalog Number	Size
PerfeCta NGS Quantification Kit - Illumina - 500 R	95154-500	500 x 20 μ l rxns
PerfeCta NGS Quantification Kit - Illumina, ROX - 500 R	95155-500	500 x 20 μ l rxns
PerfeCta NGS Quantification Kit - Illumina, Low ROX - 500 R	95156-500	500 x 20 μ l rxns

INSTRUMENT COMPATIBILITY

High Rox

- Applied Biosystems 5700
- Applied Biosystems 7000
- Applied Biosystems 7300
- Applied Biosystems 7700
- Applied Biosystems 7900
- Applied Biosystems 7900HT
- Applied Biosystems 7900HT Fast
- Applied Biosystems StepOne™
- Applied Biosystems StepOnePlus™

No Rox

- Q
- QIAGEN® Rotor-Gene® Q
- Bio-Rad CFX
- Other
- Roche Lightcycler

Low Rox

- Applied Biosystems 7500
- Applied Biosystems 7500 Fast
- Stratagene Mx3000P®
- Stratagene Mx3005P™
- Stratagene Mx4000™
- Applied Biosystems ViiA 7
- Applied Biosystems QuantStudio™
- Agilent AriaMx
- Douglas Scientific IntelliQube®

Bio-Rad iCycler iQ systems

- Bio-Rad iCycler iQ™
- Bio-Rad MyiQ™
- Bio-Rad iQ™ 5

DISCLAIMER / TRADEMARKS

References

1. Schuster, D., G. Buchman, and A. Rashtchian. (1992) A simple and efficient method for amplification of cDNA ends using 5' RACE. *Focus* 14: 46-52.
2. Buchman, G.W. and A. Rashtchian. (1992) PCR amplification of nucleic acid sequences using the 3' RACE system and direct cloning of the amplified products. *Focus* 14: 2-5.
3. Lee, E.H., Sitaraman, K., Schuster, D. and Rashtchian, A. (1997) A highly sensitive method for one-step amplification of RNA by polymerase chain reaction. *Focus* 19: 39-42
4. Westfall, B., Sitaraman, K., Solus, J., Hughes, J., and Rashtchian, A. (1997) Improved specificity and yield with Platinum Taq DNA Polymerase. *Focus* 19: 46-47.
5. Schuster, D.M., Darfler, M., Lee, J.E., and Rashtchian, A. (1998) Improved sensitivity and specificity of RT-PCR. *Focus* 20: 33-34.
6. Schwabe, W., Lee, J.E., Xu, R.H., Sitaraman, K., Smith, M., Potter, R.J., Rosenthal, K., Rashtchian, A., and Gerard, G.F. (1998) ThermoScript™ RT, A new avian reverse transcriptase for high-temperature cDNA synthesis to improve RT-PCR. *Focus* 20: 30-33
7. Westfall, B., Sitaraman, K., Lee, J., Borman, J. and Rashtchian, A. (1999) Platinum Pfx DNA Polymerase for high fidelity PCR. *Focus* 21: 46
8. Thiel, V., Rashtchian, A., Herold, J., Schuster, D.M., Guan, N., and Siddell, S.G. (1997) Effective amplification of 20-kb DNA by reverse transcription PCR. *Anal Biochem.* 252(1):62-70.
9. Xu, R.H., Schuster, D.M., Lee, J.E., Smith, M., Potter, J., Dhariwal, G., Rosenthal, K., Nathan, M., Gerard, G.F., and Rashtchian, A. (2000) One-step analysis and quantification of RNA by RT-PCR using high-temperature reverse transcription. *Focus* 22: 3-5.
10. Borman, J., Schuster, D., Li, W., Jessee, J., and Rashtchian, A. (2000) PCR from problematic template. *Focus* 22: 10-11

Trademarks

Quantabio, PerfeCra®, FastMix®, ToughMix®, qScript™, GelTrack™, and Ultraplex® are trademarks of QIAGEN Beverly, Inc..

Applied Biosystems, AmpliTaq, Gold, Gold 360, GTXpress, Path-ID™, SuperScript®, Platinum®, SYBR®, ThermoScript™, OneStep™ ViiA7, ROX™, and GreenER™ are trademarks of Life Technology Corporation.

iScript™, CFX96™, CFX384™, iQ™, MyiQ™, iCycler®, iCycler, iQ®, Chromo4™, MiniOpticon™ and Opticon® are trademarks of Bio-Rad Laboratories.

TaqMan® and LightCycler® are registered trademarks of Roche Molecular Systems, Inc.

Phase Lock Gel®, HotMaster®, Rotor-Gene® QuantiFast, and QuantiTect® are trademarks of Qiagen.

Mastercycler® is a trademark of Eppendorf.

Stratagene, MX3000P®, MX3005P™ and MX4000® are trademarks of Agilent Corporation. OneStep™ is a trademark of Applied Biosystems.

Smartcycler® is a trademark of Cepheid.

PikoReal is a trademark of Thermo Fisher Scientific.

Eco is a trademark of Illumina.

Quantabio is licensed for qPCR see www.quantabio.com for details.

© 2020 QIAGEN Beverly Inc. 100 Cummings Center Suite 407J Beverly, MA 01915

Quantabio brand products are manufactured by QIAGEN, Beverly Inc.

Version 4.0, 01/2020

Quantabio products are intended for molecular biology applications. The products are not intended for the diagnosis, prevention or treatment of a disease.

MK-PC-0002 REV 01 Catalog Quantabio 0120

Summary of Limited Label License Statements for Quantabio's PCR, qPCR, and cDNA synthesis kits

Limited Label Licenses

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this manual and for use with components contained in the kit only. QIAGEN Beverly, Inc. grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this manual, and additional protocols available at www.quantabio.com. Some of these additional protocols have been provided by Quantabio product users. These protocols have not been thoroughly tested or optimized by QIAGEN Beverly, Inc.. QIAGEN Beverly, Inc. neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN Beverly, Inc. makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN Beverly, Inc. specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN Beverly, Inc. may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

Limited Use Label License (qScript cDNA SuperMix)

This product is covered by US patent 7,470,515, US patent 7,638,612 and other patents pending in the United States and Europe. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer. The buyer is not authorized to sell or otherwise transfer this product, any of its components to a third party. The purchase of this product does not authorize the purchaser to use the product or any of its components for manufacture of commercial product. For information on obtaining a license to this product for purposes other than research, contact Licensing Department, QIAGEN Beverly Inc. 100 Cummings Center Suite 407J Beverly, MA 01915.

Limited Label License (AccuStart II, GelTrack, ToughMix and FastMix II products).

This product is licensed under US 7,972,828 and corresponding US and foreign patents and patent applications.



A faster, smaller, better way to qPCR



Ultra-Fast Data Acquisition 35 cycles in 25 minutes



Unrivaled Performance Detect two-fold differences in expression levels reproducibly between samples, runs and instruments



Portable & Compact At 4.5 lbs & speaker-sized, transport without ever calibrating



Scalable & Wireless Wirelessly connect & operate up to 10 instruments from one PC (48 samples/instrument)



Magnetic Induction Technology Temperature accuracy of $\pm 0.25^{\circ}\text{C}$ combined with $\pm 0.05^{\circ}\text{C}$ well-to-well uniformity eliminates variability vs block-based cyclers

Description

Q uses a patented magnetic induction technology to rapidly heat samples coupled with fan forced air for cooling to acquire data in only 25 minutes. The robust optical system acquires all channels simultaneously and allows for running the fastest multiplexed assays.

Q's miniature speaker-size and 4.5 pound weight make it the most portable and versatile qPCR cyclers on the market without ever needing to calibrate. Q also provides scalability as each instrument can process up to 48 samples per run and up to 10 Q's can be connected to a single computer wirelessly via bluetooth enabling up to 480 samples to be processed simultaneously.

A key difference is that Q incorporates a unique spinning aluminum rotor providing superior temperature uniformity of $\pm 0.05^{\circ}\text{C}$ versus traditional block-based real time cyclers which rely on multiple peltier heating blocks that can create edge effects resulting in sample variation. Not only does the data give you superior reproducibility, repeatability but enables detection of 2-fold gene expression level differences as well as identification of difficult class IV SNP's requiring melt temperature resolutions of 0.1°C .

Who wouldn't want to take one for a spin?



25 mins

Ultra-Fast Data Acquisition

Generate high quality data, fast!

Q's speed is the fastest in the industry. Coupling speed with superior temperature uniformity means you don't sacrifice on the performance quality of your qPCR. Completing runs in as little as 25 minutes* is the new standard.

Normalized Fluorescence





Figure 1

5 point, 2x dilution series of Hepatitis B virus (HBV) cDNA template
Starting amount of $3\text{E}+06$ copies ($n = 4$ each)
Efficiency = 90% (standard curve method); $R^2 = 0.99$
Time to complete run (including melt) = 26 min


*25 minute cycle times obtained with fast cycling master mixes and short amplicon assay designs targeting cDNA



Unrivaled Performance

Detect 2-fold Differences

Normalized Fluorescence




Confidently detect small differences.

Attributed to the high thermal uniformity and reproducibility, Q is able to reliably detect 2-fold differences in gene expression levels whether using standard curves for absolute quantification or relative quantification through REST.

Figure 2

Manganese superoxide dismutase gene (MnSOD)
Eight point, 2x dilution series of human genomic DNA (n = 4 each)
Efficiency = 98% (standard curve method)
 $R^2 = 1.00$

Normalized Fluorescence



Q is also able to detect differences within a single cycle. As an example, the data shown (left) illustrates a 5 pg dilution series differentiated with 0.2 cycles between standards.

Figure 3

Five point dilution series of HBV plasmid cDNA template (n = 4 each)
5 picogram difference between standards
Efficiency = 98% (standard curve method)
 $R^2 = 0.99$

Wide Linear Dynamic Range

Q provides ultra-sensitive detection down to single digit copies of DNA. Whether using absolute quantification for viral load detection or a standard curve to determine PCR efficiency, Q performs reliably across a wide dynamic range.

Normalized Fluorescence




Figure 4

10 point, 10x dilution series of Hepatitis B virus (HBV) cDNA template
Starting amount of $3E+09$ copies (n = 3 each) over 10 logs
Efficiency = 95% (standard curve method)
 $R^2 = 0.99$





Portable & Compact

4.5 lbs

A qPCR cycler the size of a speaker!

Q is the most portable and compact cycler on the market occupying $\frac{1}{4}$ the size of current qPCR platforms. The small footprint is enabled by the unique rotary and magnetic induction technology.



Don't be fooled by it's small size. Big things are happening inside.

- Fixed optics & no moving parts
- Never needs optical alignment or calibration
- No reference dyes or cross talk compensation required
- Utilizes proprietary 0.1 ml 4-strip tubes & caps supporting volumes of 5 – 30 μ l
- High speed centrifugation ensures sample spin down
- Prevents evaporation & condensation with pre-loaded oil in tubes



up to 10 cyclers

Scalable & Wireless


Flexibility from 48 to 480 samples

Each Q cycler can process up to 48 samples per run and up to 10 cyclers can be connected to a single computer wirelessly to provide the desired scalability.

Q's advanced data analysis software makes combining a single data set from multiple runs from multiple cyclers seamlessly simple.

High reproducibility from thermal and optical performance ensures data from different runs and cyclers look like it was generated on the same instrument and the same run.

Lastly, anyone can setup and run the Q as it is plug-and-play right of the box.




Specifications



Physical

Height	5.1 in
Width	5.9 in
Length	5.9 in
Weight	4.5 lbs



Thermal Performance

Temperature Accuracy	$\pm 0.25^{\circ}\text{C}$
Temperature Uniformity	$\pm 0.05^{\circ}\text{C}$
Ramp Rates	Heating 4°C/s / Cooling 3°C/s
Temperature Input Range	40 – 99 $^{\circ}\text{C}$

Optical

Detectors	Photodiode per channel		
Excitation Sources	High power LED for each channel		
Channels	Green	Ex 465 nm	Em 510 nm
	Yellow	Ex 540 nm	Em 570 nm
	Orange	Ex 585 nm	Em 618 nm
	Red	Ex 635 nm	Em 675 nm
Acquisition Time	1 second		

Reaction Tubes

Samples per Instrument	48
Reaction Volume Range	5 – 30 μl

Operating Environment

Temperature	18 – 35 $^{\circ}\text{C}$
Relative Humidity	20 – 80%

ORDER INFO

Product Name	Quantabio Catalog Number	Size
Q 2-channel qPCR Instrument	95900-4C	1 instrument
Q Tubes & Caps (20 racks/box)	95910-20	1 box
1 year extended warranty 4-channel Q Instrument	95900-4CEXTWAR	–
HRM License for Q qPCR Instrument	95915-20	1 HRM License
Related Products		
SYBR based qPCR		
PerfeCta® SYBR Green FastMix®	95072-012	10 x 1.25 ml
PerfeCta SYBR Green FastMix	95072-05K	1 x 50 ml
Probe based qPCR		
PerfeCta qPCR FastMix II	95118-012	10 x 1.25 ml
PerfeCta qPCR FastMix II	95118-05K	1 x 50 ml
Genotyping ToughMix®	95115-05K	1 x 50 ml
PerfeCta qPCR ToughMix	95112-012	10 x 1.25 ml
PerfeCta qPCR ToughMix	95112-05K	1 x 50 ml
NGS Library Quantification		
sparQ Universal Library Quant Kit - 100	95210-100	100 rxns
sparQ Universal Library Quant Kit - 500	95210-500	500 rxns
Quantitative RT-qPCR		
qScript® One-Step SYBR Green RT-qPCR Kit	95087-200	4 x 1.25 ml
qScript XLT One-Step RT qPCR ToughMix	95132-100	1 x 1 ml
qScript XLT One-Step RT qPCR ToughMix	95132-500	5 x 1 ml
qScript One-Step RT-qPCR Kit	95057-200	4 x 1.25 ml
UltraPlex One-Step ToughMix	95166-100	1 x 500 µl
UltraPlex One-Step ToughMix	95166-500	5 x 500 µl

Q Cyclers does not require the use of reference dyes.

CONTACT



100 Cummings Center
Beverly, MA 01915
USA

For a comprehensive listing of our distribution partners please visit our website.

www.quantabio.com

